



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com


Transparent Female Catheter Model

Item No. 1017242 [W43028F]
Weight 9.53kg
Dimensions 52.7 x 24.1 x 41 cm
MPN: M160

[Read More](#)

SKU:

Categories: Catheterization



Product Description

Makes it possible to practice catheterization while observing what is happening internally. Allows you to view the indwelling condition of a balloon catheter, flow of urine, and the relative position of the pelvis and bladder through a transparent hypogastrium. External genitalia is made of a special soft material,

making this model extremely realistic. Comes with irrigator, irrigator stand, and case

Why Choose the Transparent Female Catheter Model for Your Nursing Skills Lab in India?

The Transparent Female Catheter Model (Item No. 1017242 [W43028F], 9.53 kg, 52.7 × 24.1 × 41 cm, MPN: M160) lets students practice catheterization while directly observing what's happening internally — a genuinely distinct teaching approach from fluid-feedback-only trainers.

Through a transparent hypogastrium, students can watch the indwelling condition of the balloon catheter, the actual flow of urine, and the relative position of the pelvis and bladder as the procedure unfolds. This visual insight helps students build an accurate mental model of the internal anatomy involved, not just the external procedural steps. External genitalia are made of a special soft material for a realistic feel during insertion.

Includes: Irrigator, irrigator stand, carrying case.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of leading simulation training equipment in India, with over 30 years of experience supplying medical simulators and clinical skills training equipment to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: What can I observe through the transparent hypogastrium? A: The indwelling position of the balloon catheter, urine flow, and the relative position of the pelvis and bladder in real time.

Q: How do I order from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.