



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



Life/form® Single Intraosseous Infusion Leg

Item No. LF03622U
Weight 0.8kg
Dimensions 30 x 25 x 10 cm
Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories: Injections and Punctures

Product Description

This simulator is a single leg with the ability to be placed on the Life/form® Infant CRiSis™ Manikin. The leg comes complete with 10 replaceable bones, two sets of replacement skins, and a fluid administration kit enabling you to actually visualize a flashback and infuse fluid into the leg. Intraosseous infusion needle included. Size: 12" x 10" x 4".

Why Choose the Life/form® Single Intraosseous Infusion Leg for Your Nursing Skills Lab in India?

The Life/form® Single Intraosseous Infusion Leg (Item No. LF03622U, Weight: 0.8kg, Dimensions: 30 x 25 x 10 cm, Make: Nasco Healthcare, USA) is built for practicing intraosseous (IO) access in the leg — the standard first-line site for pediatric IO infusion in real emergencies, distinct from the sternal

approach used in adult IO training.

This single leg simulator is designed to be placed on the Life/form® Infant CRiSis™ Manikin, letting instructors integrate IO access practice into broader infant emergency scenario training.

What's included:

- 10 replaceable bones, allowing repeated puncture practice without needing to replace the whole unit each time
- Two sets of replacement skins
- A fluid administration kit that lets students visualize true flashback and infuse fluid into the leg
- Intraosseous infusion needle

Specifications:

- Size: 12" x 10" x 4"

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators and clinical skills training equipment to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: Which manikin is this leg designed to attach to? A: The Life/form® Infant CRiSis™ Manikin.

Q: How many times can the leg be punctured before needing new parts? A: It includes 10 replaceable bones and two sets of replacement skins, allowing extended reuse before full replacement is needed.

Q: Does this simulator show a realistic flashback? A: Yes — the included fluid administration kit lets students visualize flashback and infuse fluid into the leg upon successful IO access.

Q: How do I order from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.