



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com

Basic GERi™ Nursing Skills Elderly-Care Manikin

Item No. 1005620 [W44077]
Weight 10.785kg
Dimensions 94 x 52 x 28 cm
MPN: LF04040U

[Read More](#)

SKU:

Categories: Intramuscular (I.m.) and Intradermal



Product Description

Why Choose the Basic GERi™ Nursing Skills Elderly-Care Manikin for Your Simulation Lab in India?

The Basic GERi™ Nursing Skills Elderly-Care Manikin (Item No. 1005620 / W44077, MPN: LF04040U)

by Nasco Healthcare, USA is a full-size, lightweight elderly-care simulation manikin weighing 10.785 kg and measuring 94 × 52 × 28 cm. Designed with realistic aged skin — including wrinkles, skin folds, and age-appropriate anatomical proportions — GERi™ provides nursing students and healthcare trainees with a lifelike platform for practising the full breadth of basic elderly patient care skills. For NMC-mandated clinical skills labs, geriatric nursing training programmes, hospital in-service education teams, and nursing colleges across India seeking a dedicated geriatric manikin at the Basic tier, GERi™ is the entry point into the Life/form® GERi™ range.

India's ageing population places increasing demands on nurses and caregivers who must be competent in elderly-specific care. Training on a manikin that accurately replicates the physical characteristics of an older adult — including reduced skin elasticity, realistic joint range of motion, and age-appropriate anatomy — prepares trainees for the realities of caring for elderly patients in ways that a standard adult manikin cannot.

The Basic GERi™ supports a wide range of training scenarios. For visual inspection training, the manikin presents a normal and cancerous mole, a Stage 1 sacral pressure ulcer, dilated and constricted pupils, and reddened skin folds — giving trainees practice in identifying clinically significant skin and ocular findings. For patient care simulation, the manikin supports bandaging and wound dressing, bed baths, clothing changes, denture placement and removal (upper and lower), ear canal irrigation, otic drop administration and hearing aid placement, eye irrigation, finger and toe manipulation with bandaging, hair care (washing and combing), intramuscular injection at arm, thigh, and buttock sites, oral and nasal hygiene including lavage, gavage and suctioning, ostomy care covering ileostomy and colostomy tissue maintenance and appliance application, patient positioning, patient transfer techniques, pericare, and tracheostomy care with lavage and suctioning. The manikin presents as female by default and converts to male with removal of the wig and attachment of the included male genitalia. Joints are non-pinching and move like a real person, supporting realistic patient positioning and transfer practice. The Basic GERi™ is ideal for OBRA-required training. Carry case is not included.

For institutions requiring advanced capabilities, GERi™ is also available in Complete (adds urinary catheterisation, nasogastric tube placement, gastrostomy, enema, pap smear, and prostate exam), Advanced (adds IV insertion and blood pressure monitoring), and Auscultation (adds heart and lung sound recognition) configurations.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators to 1,500+ institutions across 28 states. We provide pan-India delivery, on-site installation, faculty orientation, and dedicated after-sales support.

Q: Is this manikin male or female? A: The Basic GERi™ presents as female by default. It converts to male with removal of the wig and attachment of the male genitalia included in the package.

Q: What skin findings does GERi™ display for visual inspection training? A: The manikin features a normal mole, a cancerous mole, a Stage 1 sacral pressure ulcer, dilated and constricted pupils, and reddened skin folds.

Q: What is the difference between Basic, Complete, Advanced, and Auscultation GERi™? A: The Basic version covers the core nursing skills listed above. The Complete adds urinary catheterisation, nasogastric tube placement, gastrostomy, enema administration, pap smear, and prostate exam. The Advanced adds IV insertion and blood pressure monitoring. The Auscultation model adds heart and lung sound recognition. SEM Trainers can supply all four configurations — contact us to discuss which tier suits your curriculum.

Q: Is a carry case included? A: No. A carry case is not included with the Basic GERi™

GERi™ has an elderly appearance with skin wrinkles and folds. Complete lifelike range of motion, realistic patient positioning, and non-pinching joints - moves like a person. • Full-size adult manikin - measures 58". • Lightweight - approximately 28 lbs. • Overall female appearance • Converts to male with removal of wig and attachment of male genitals
Visual Inspection: • Normal and Cancerous Mole • Stage 1 Sacral Ulcer • Dilated and Constricted pupils • Reddened Skin Folds
Patient Care Simulation: • Bandaging and Wound Dressing • Bed Baths • Clothing Changes • Denture Placement and Removal, Upper and Lower • Ear Canal Irrigation, Otic Drops, and Hearing Aid Placement • Eye Irrigation • Finger and Toe Manipulation with Bandaging - Flexible and Individually Molded • Hair Care Washing and Combing • Intramuscular Injection Sites: Arm, Thigh, Buttock • Oral and Nasal Hygiene: Lavage, Gavage, and Suctioning • Ostomy Care: Ileostomy and Colostomy Tissue • Maintenance and Appliance Application, Lavage, and Suctioning • Patient Positioning • Patient Transfer Techniques • Pericare • Tracheostomy Care, Lavage and Suctioning Ideal for OBRA Required Training. Carry case not included.

	Basic: KERi™ 1020187 GERi™ 1005620	Complete: KERi™ 1013742 GERi™ 1005597	Advanced: KERi™ 1017567 GERi™ 1005606	Auscultation: KERi™ 1020147 GERi™ 1020146
Skin inspection and care	Yes	Yes	Yes	Yes
Bandaging	Yes	Yes	Yes	Yes
Washing	Yes	Yes	Yes	Yes
Denture Placement	Yes	Yes	Yes	Yes
Hearing aid placement	Yes	Yes	Yes	Yes
Ear and eye irrigation	Yes	Yes	Yes	Yes
Ostomy and tracheostomy care	Yes	Yes	Yes	Yes
Positioning and transfer techniques	Yes	Yes	Yes	Yes
Pericare	Yes	Yes	Yes	Yes
Intramuscular injections	Yes	Yes	Yes	Yes

Gastrostomy	-	Yes	Yes	Yes
Pap smears and douching	-	Yes	Yes	Yes
Enema administration	-	Yes	Yes	Yes
Prostate exam	-	Yes	Yes	Yes
Urinary catheterisation	-	Yes	Yes	Yes
Nasogastric tube placement	-	Yes	Yes	Yes
IV insertion	-	-	Yes	
Blood pressure monitoring	-	-	Yes	
Recognition of heart sounds	-	-	-	Yes
Recognition of lung sounds	-	-	-	Yes