



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Prostate Examination Trainer

Item No. 1005163 [W19344]
Weight 0.11kg
MPN: BPH100

[Read More](#)

SKU:

Categories: Male Examination

Product Description

The prostate examination trainer is supplied as a four sided box housing 4 prostates: Normal Benign, Enlarged Benign, Early Malignant, Advanced Malignant. Features of the prostate examination trainer include:

- Realistic tactile feeling
- No prostates need to be replaced
- Each is identified via a single, gloved finger access port
- Easy to handle

The soft but firm touch of these prostates has been extensively tried and tested by specialist urologists. The prostate examination trainer is great for all types of prostate education.

Why Choose the Prostate Examination Trainer (BPH100) for Your Clinical Skills Programme in India?

The Prostate Examination Trainer (Item No. 1005163 / W19344, MPN: BPH100) is a compact, zero-maintenance prostate palpation familiarisation trainer from Pharmabotics, UK / 3B Scientific — designed for medical students, urology trainees, general practice residents, and family medicine students in India who need to develop the tactile vocabulary of prostate palpation before they encounter the findings on a real patient or move to a full DRE technique simulator.

Weighing just 0.11 kg, the BPH100 is the most portable and accessible prostate palpation reference tool available — a four-sided box housing four individual prostates, each accessible through its own dedicated gloved finger access port. No body form, no lubricant, no replacement parts, no setup. The trainee inserts one gloved finger through the access port and palpates the prostate — exactly as described, nothing more — building the discriminatory touch and descriptive language that makes a DRE clinically useful.

Four Prostate Presentations — Core Clinical Spectrum:

Each of the four sides of the box presents one prostate for palpation:

- **Normal Benign** — baseline reference for normal size, shape, and consistency
- **Enlarged Benign (BPH)** — symmetric enlargement without nodularity
- **Early Malignant** — palpable asymmetry or nodularity on an otherwise firm gland
- **Advanced Malignant** — hard, irregular, fixed presentation characteristic of advanced prostate carcinoma

The soft but firm tactile properties of each prostate have been developed and validated by specialist urologists, ensuring the palpatory experience is clinically credible and representative of real findings.

No Replacements Required: Unlike simulators with consumable inserts, the BPH100 prostates are durable and require no replacement — making this a zero-ongoing-cost training tool suitable for high-volume student throughput in skills labs, lectures, and OSCE preparation sessions.

Where the BPH100 Fits in the Prostate Training Progression:

The BPH100 is the entry point of the prostate training pathway — ideal for the first encounter with prostate palpation feel and for building pathology description skills. For institutions requiring full DRE technique training including approach, insertion angle, and patient positioning, SEM Trainers also supplies the Advanced Patient Care Male Prostate Simulator (1020127, Gaumard Scientific) and the Prostate Examination Simulator M53B (Kyoto Kagaku, 10 cases). See the comparison table below or

contact SEM Trainers to discuss which product fits your programme level.

Why Buy from SEM Trainers? SEM Trainers & Systems has over 30 years of experience supplying medical simulators to 1,500+ institutions across 28 states in India. We provide pan-India delivery, on-site installation, faculty orientation, and dedicated after-sales support.

Q: How does the trainee access each prostate? A: Each of the four sides of the box has a dedicated gloved finger access port — one per prostate. The trainee inserts one gloved finger through the port and palpates the prostate directly. No lubricant is required, though a glove is recommended for hygiene.

Q: Do the prostates ever need to be replaced? A: No. The four prostates are durable and integrated into the box — no consumable replacements are required. This makes the BPH100 a zero-ongoing-cost training tool for high-volume student use.

Q: What is the difference between this trainer and the Prostate Examination Simulator M53B? A: The BPH100 is a compact palpation familiarisation tool — four prostates in a portable box, no body form, no insertion anatomy. The M53B (Kyoto Kagaku) is a premium departmental DRE training suite with 10 prostate cases on a rotating internal table, a full male lumbar torso body form, and three patient position options. The BPH100 trains palpation feel and description; the M53B trains the complete DRE procedure. Both are available from SEM Trainers.

Q: Is this trainer suitable for OSCE station use? A: Yes — its compact size and zero-setup design make it ideal for OSCE prostate palpation description stations, peer-to-peer practice, and lecture demonstrations.

Q: What is the difference between this trainer and the Advanced Patient Care Male Prostate Simulator (1020127)? A: The BPH100 trains prostate palpation feel only — no insertion anatomy. The 1020127 (Gaumard Scientific) includes a realistic pelvic body form with anatomically correct anal opening and four card-mounted prostate modules, training the full DRE technique including approach and insertion. The BPH100 is the affordable starting point; the 1020127 is the step-up to full procedural DRE training.