



**SEM**

## SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

**Mobile :** +91 88495 63724

**Mobile 1:** +91 98791 03905

**Email:** sem@semtrainers.com

**Website:** www.semtrainers.com



## Stomach Bag Replacement for Keri / Geri

Item No. 1019751  
Weight 0.01kg  
MPN: LF04098U

[Read More](#)

**SKU:**

**Categories:**Adult Patient Care

## Product Description

Component replacement part for KERi™/GERi™ Nursing Skills Manikin.

**Stomach Bag Replacement for KERi™/GERi™ Nursing Skills Manikin (1019751 / LF04098U) — Genuine OEM Replacement Consumable**

The Stomach Bag Replacement for KERi™/GERi™ (Item No. 1019751, MPN: LF04098U) is a genuine Nasco Healthcare Life/form® replacement gastric reservoir bag for the KERi™ and GERi™ Nursing Skills Manikin range. It is the internal stomach bag housed inside the manikin torso — the fluid reservoir that enables nasogastric (NG) tube feeding, orogastric (OG) tube feeding, gastric lavage, and gastric suction training. It is the correct reorder item when the internal stomach bag on your KERi™ or GERi™ manikin has degraded, perforated, or become unusable through repeated use.

**Compatible with:** KERi™ and GERi™ Nursing Skills Manikins — all tiers (Basic, Complete, Advanced) that include NG tube, OG tube, or gastric suction training capability. Confirm compatibility with your specific manikin model before ordering.

**Why Buy from SEM Trainers?** SEM Trainers & Systems has over 30 years of experience supplying medical simulators and genuine spare parts to 1,500+ institutions across 28 states in India. We provide pan-India delivery and dedicated after-sales support.

### Frequently Asked Questions

**Q: What is the difference between the Stomach Bag Replacement (1019751) and the Stoma Bag (LF04099U)?** A: The Stomach Bag Replacement (1019751 / LF04098U) is the internal gastric reservoir inside the manikin torso — used for NG tube, OG tube, gastric lavage, and gastric suction training. The Stoma Bag (LF04099U) is the external ostomy pouch that attaches over the stoma opening — used for colostomy and ileostomy care training. They are completely different components. If you are unsure which part you need, contact SEM Trainers with your manikin's item number.

**Q: Which manikins is this compatible with?** A: KERi™ and GERi™ Nursing Skills Manikins from Nasco Healthcare — all tiers that include NG/OG tube and gastric suction training. Confirm compatibility with your specific manikin model before ordering.

**Q: Is this a genuine OEM part?** A: Yes. Item No. 1019751 (MPN: LF04098U) is a genuine Nasco Healthcare Life/form® replacement consumable.

**Q: When does the stomach bag need replacing?** A: When the internal bag shows leakage, perforation, or significant degradation that prevents it from holding simulated fluid reliably during NG/OG tube and gastric suction training sessions.