



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com

Advanced GERi™ Manikin

Item No. 1005606 [W44046]
Weight 27kg
Dimensions 95.5 x 52 x 28 cm
MPN: LF04030U

[Read More](#)

SKU:

Categories: Geriatric Patient Care



Product Description

GERi™ has an elderly appearance with skin wrinkles and folds. Complete lifelike range of motion, realistic patient positioning, and non-pinching joints - moves like a person. • Lightweight - approximately 28 lbs. • Full-size adult manikin - measures 58" • Overall female appearance • Converts to male with removal of wig and attachment of male genitals Visual Inspection: • Normal and Cancerous Mole • Stage 1 Sacral Ulcer • Dilated and constricted pupils • Reddened Skin Folds Patient Care Simulation: • Bandaging and Wound Dressing • Bed Baths • Clothing Changes • Denture Placement and Removal, Upper and Lower • Ear Canal Irrigation, Otic Drops, and Hearing Aid Placement • Eye Irrigation • Finger and Toe Manipulation with Bandaging - Flexible and Individually Molded • Hair Care Washing and Combing • Intramuscular Injection Sites: Arm, Thigh, Buttock • Oral and Nasal Hygiene: Lavage, Gavage, and Suctioning • Ostomy Care: Ileostomy and Colostomy Tissue • Maintenance and Appliance Application, Lavage, and Suctioning • Patient Positioning • Patient Transfer Techniques • Pericare •

Tracheostomy Care, Lavage and Suctioning • Gastrostomy Procedures, Lavage and Gavage • Enema Administration, Female Only • Urinary Catheterization, Female and Uncircumcised Male • Pap Smears and Douching • Prostate Exam, Stage B • Nasogastric Tube Placement • Internal Fluid Reservoirs The Advanced GERi™ Manikin simulates over 35 nursing and medical procedures, and also includes: • Right IV Training Arm with a shoulder intramuscular injection site featuring replaceable skin and veins that roll when palpated; realistic flashback confirms proper needle placement • Left Blood Pressure Training Arm that reproduces all five Korotkoff sounds and allows the instructor to vary systolic and diastolic levels, pulse rate, volume, and auscultatory gap • Standard arms with right and left shoulder intramuscular injection sites are also included.

Advanced GERi™ Manikin (1005606 [W44046] / LF04030U)

GERi™ has an elderly appearance with skin wrinkles and folds, offering complete lifelike range of motion for realistic patient positioning and non-pinching joints that truly move like a person.

A true full-size, lightweight adult manikin. GERi™ measures 58" and weighs approximately 28 lbs, with an overall female appearance that easily converts to male with removal of the wig and attachment of male genitalia.

Visual inspection findings built right in: normal and cancerous mole, a stage 1 sacral ulcer, dilated and constricted pupils, and reddened skin folds.

Simulates over 35 nursing and medical procedures, including bandaging and wound dressing, bed baths, clothing changes, denture placement and removal, ear canal irrigation with otic drops and hearing aid placement, eye irrigation, finger and toe manipulation with bandaging, hair care, intramuscular injection sites (arm, thigh, buttock), oral and nasal hygiene (lavage, gavage, suctioning), ostomy care with appliance maintenance, patient positioning and transfer techniques, pericare, tracheostomy care, gastrostomy procedures, enema administration (female), urinary catheterization (female and uncircumcised male), pap smears and douching, prostate exam (Stage B), nasogastric tube placement, and internal fluid reservoirs.

Advanced tier adds dedicated IV and blood pressure training arms. A right IV training arm features a shoulder intramuscular injection site with replaceable skin and veins that roll when palpated, plus realistic flashback to confirm proper needle placement. A left blood pressure training arm reproduces all five Korotkoff sounds, letting instructors vary systolic and diastolic levels, pulse rate, volume, and the auscultatory gap. Standard arms with right and left shoulder intramuscular injection sites are also included.

Ideal for OBRA-required training. Carry case not included.

Part of a tiered GERi™/KERi™ family. This Advanced tier builds on Basic and Complete skin care, ostomy, and gastrostomy/catheterization skills by adding dedicated IV insertion and blood pressure

monitoring training arms, while the Auscultation tier focuses on heart and lung sound recognition.

Why Buy from SEM Trainers? SEM Trainers & Systems has over 30 years of experience supplying medical education resources to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: What does the Advanced tier add over Complete? A: A dedicated right IV training arm with realistic flashback confirmation, and a left blood pressure training arm reproducing all five Korotkoff sounds.

Q: Is a carry case included? A: No — the carry case is not included and must be purchased separately.

Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.

Ideal for OBRA Required Training. Carry case not included.

	Basic: KERi™ 1020187 GERi™ 1005620	Complete: KERi™ 1013742 GERi™ 1005597	Advanced: KERi™ 1017567 GERi™ 1005606	Auscultation: KERi™ 1020147 GERi™ 1020146
Skin inspection and care	Yes	Yes	Yes	Yes
Bandaging	Yes	Yes	Yes	Yes
Washing	Yes	Yes	Yes	Yes
Denture Placement	Yes	Yes	Yes	Yes
Hearing aid placement	Yes	Yes	Yes	Yes
Ear and eye irrigation	Yes	Yes	Yes	Yes
Ostomy and tracheostomy care	Yes	Yes	Yes	Yes
Positioning and transfer techniques	Yes	Yes	Yes	Yes
Pericare	Yes	Yes	Yes	Yes
Intramuscular injections	Yes	Yes	Yes	Yes
Gastrostomy	-	Yes	Yes	Yes
Pap smears and douching	-	Yes	Yes	Yes
Enema administration	-	Yes	Yes	Yes
Prostate exam	-	Yes	Yes	Yes
Urinary catheterisation	-	Yes	Yes	Yes
Nasogastric tube placement	-	Yes	Yes	Yes
IV insertion	-	-	Yes	
Blood pressure monitoring	-	-	Yes	
Recognition of heart sounds	-	-	-	Yes

Recognition of lung sounds	-	-	-	Yes
----------------------------	---	---	---	-----