



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Standard Ready-or-Not Tot® - White Male

Item No. 1017931 [W44214]
Weight 2.988kg
Dimensions 56.5 x 29.6 x 16.5 cm
MPN: WA18306

[Read More](#)

SKU:

Categories: Neonatal Patient Care



Product Description

The original manikin that simulates a baby's varying needs! Take the glamour out of teen pregnancy and demonstrate the realities of the real world to your students. Characteristics of Standard Ready-or-Not Tot®: • Offers three different tending programs typical of most babies, to simulate normal baby-tending needs along with fussy periods in each program. • Provides programmed tending events

that are unpredictable to the student, but easy for the teacher to monitor, unlike other simulators. • Cries, coos, and burps in response to appropriate care. Programs are 48 hours in length and run continuously. • Includes a “demo” option. • Demands student’s attention for periods of 5 minutes to 30 minutes for each tending occurrence. There are 25-27 tending occurrences within each 48-hour time frame (for total of 7 hours of care). • Provides teacher with complete control of the child care experience. • Includes a “panic” key for quieting difficulties without discontinuing the entire simulation. • Alerts teacher to a student’s abuse, panic, or tampering with the control box. • Newborn-size doll represents a baby of a teen mom - 5 lbs. 8 oz., 20" long. • Anatomically correct. You never have to worry if students neglected the manikin. Ready-or-Not Tot® uses set programs so the teacher always knows when the manikin required care. Student response to the manikin’s needs are easily checked using the template provided. Comes complete with one set of student keys, one set of teacher keys, reproducible parental consent form/permission slip, reproducible student response sheet, teacher correction template, diaper, 9V battery, and instructions.

Standard Ready-or-Not Tot® – White Male (1017931 [W44214])

The Standard Ready-or-Not Tot® – White Male (Item No. 1017931 [W44214], 2.988 kg, 56.5×29.6×16.5 cm, MPN WA18306) is the original manikin that simulates a baby's varying care needs, giving students a realistic look at the demands of infant care — a valuable tool for teen pregnancy prevention education.

Three realistic tending programs. The Tot® offers three different tending programs typical of most babies, simulating normal baby-tending needs along with fussy periods within each program.

Unpredictable for students, transparent for teachers. Programmed tending events are unpredictable to the student but easy for the teacher to monitor — a key advantage over other simulators. The Tot® cries, coos, and burps in response to appropriate care, running continuously on 48-hour programs, with a "demo" option also included.

Realistic demands built into every session. Each tending occurrence demands the student's attention for 5 to 30 minutes, with 25–27 occurrences within each 48-hour period, totaling roughly 7 hours of hands-on care.

Full teacher oversight and safety features. Teachers retain complete control of the child care experience, with a "panic" key to quiet difficulties without discontinuing the simulation, and automatic alerts for suspected abuse, panic, or control-box tampering.

Realistic newborn size, representing a baby of a teen mother at 5 lbs. 8 oz. and 20" long, and anatomically correct.

No guesswork on student compliance. Since Ready-or-Not Tot® runs on set programs, teachers

always know exactly when the manikin required care, and student response is easily checked against the provided template.

Comes complete with: one set of student keys, one set of teacher keys, a reproducible parental consent form/permission slip, a reproducible student response sheet, a teacher correction template, a diaper, a 9V battery, and instructions.

Why Buy from SEM Trainers? SEM Trainers & Systems has over 30 years of experience supplying medical education resources to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: How long does each program run? A: 48 hours, continuously, with 25–27 tending occurrences totaling roughly 7 hours of hands-on care.

Q: What happens if a student mishandles the manikin? A: The system automatically alerts the teacher to suspected abuse, panic, or control-box tampering.

Q: What's included? A: Student and teacher key sets, a parental consent form, a student response sheet, a teacher correction template, a diaper, a 9V battery, and instructions.

Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.