



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Tactical Hemorrhage Control Trainer - THCT

Item No. 1021644
Weight 82.5 kg
Dimensions 192 x 70 x 20 cm
Brand Operative Experience Inc.
MPN: THCT Standard

[Read More](#)

SKU:

Categories:TCCC Training Manikins



Product Description

The Tactical Hemorrhage Control Manikin is a full-sized, remotely-activated simulator designed for point-of-injury, tactical medicine training for law enforcement and first responders. Ultra-realistic and fully mobile, it is ideally suited for training rapid assessment and treatment of trauma injuries in disaster drill and active shooter simulation scenarios.

Full-body Manikin for TCCC and TECC Training

- Anatomically-accurate soft tissue
- Ruggedized and water-resistant
- Remotely-activated pulsatile bleeding
- Multiple injuries involving gunshot wounds, stab wounds and crushing injuries
- Amputation on left leg above knee for tourniquet application

Rapid Assessment for Tactical Casualty Care

- Airway
- Oral and nasal
- Bleeding

Airway Management Training

- Oropharyngeal and Nasopharyngeal airway
- Chest Seal application
- Practice the identification of the proper landmarks for needle decompression (mid-clavicular line / mid-axillary line)

Hemorrhage Control and Tourniquet Application Training

- Gunshot wound to face with exit wound
- Gunshot wound to upper left arm with exit wound for wound packing and tourniquet application
- Superficial gunshot wounds to the chest and stomach
- Junctional and standard tourniquet application
- Practice wound packing and bandaging technic

Watch Mick Navin, CEO of Operative Experience, Inc. on TV using a [TCCC Manikin](#) to explain

hemorrhage control and tourniquet application techniques for "Stop the Bleed" day: