



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com

Episiotomy and Suturing Trainer

Item No. 1019639 [P95]
Weight 0.553kg
Dimensions 20.5 x 13.5 x 6.5 cm
Brand 3B Scientific

[Read More](#)

SKU:

Categories:Obstetrics



Product Description

Obstetrics simulation to improve clinical skills: The episiotomy and suture trainer provides realistic help with practicing taking care of and closing birth injuries after vaginal births. The medical personnel (midwives and doctors) should be continuously trained in diagnosing and taking care of birth injuries. With this trainer, taking care of the following birth injuries can be practiced and demonstrated:

- Episiotomy, for example:
 - Median episiotomy
 - Mediolateral episiotomy
 - Lateral episiotomy
- Vaginal tears
- Labial tears

To make caring for and closing these birth wounds as realistic as possible, the following muscle structures can be palpated and sutured in the perineal area:

- M. sphincter ani externus (outer anal sphincter muscle)
- M. transversus perinei superficialis (superficial transverse perineal muscle)
- M. bulbospongiosus (bulbospongiosus muscle)

Muscle and skin layers can be sutured separately. Each cut can be sutured several times; by repeating the procedure, the midwifery and obstetrics students learn about professional wound management. The suture technique can be recorded and assessed by the trainers so that students can improve, after which the students can practice again once the sutures have been removed. The material used for the simulation of the tissue feels realistic, and it is easy to keep clean. The complete trainer does not take up much room when stowed away. The adjustable holder for the simulated tissue stands firmly on the table without slipping. The insert can be taken out of the holder and then worked on while it is lying on the table. This value-for-money tissue insert can be bought individually or in packs of 5 in our shop.

Delivery contents:

- Silicone insert
- Base
- Stand

Latex-free, silicone-based material

Episiotomy and Suturing Trainer (1019639 [P95])

The episiotomy and suture trainer (Item No. 1019639 [P95], 0.553 kg, 20.5×13.5×6.5 cm) by 3B Scientific provides realistic obstetrics simulation for practicing the care and closing of birth injuries after

vaginal births. Medical personnel — midwives and doctors — should be continuously trained in diagnosing and caring for birth injuries.

A comprehensive range of birth injuries can be practiced and demonstrated, including episiotomy (median, mediolateral, and lateral), vaginal tears, and labial tears.

Realistic muscle-level practice for authentic wound care. To make caring for and closing these birth wounds as realistic as possible, the following muscle structures can be palpated and sutured in the perineal area: M. sphincter ani externus (outer anal sphincter muscle), M. transversus perinei superficialis (superficial transverse perineal muscle), and M. bulbospongiosus (bulbospongiosus muscle). Muscle and skin layers can be sutured separately.

Built for repeated, assessable practice. Each cut can be sutured several times, letting midwifery and obstetrics students learn professional wound management through repetition. Suture technique can be recorded and assessed by trainers so students can improve, then practice again once sutures are removed.

Practical, space-saving design. The tissue material feels realistic and is easy to keep clean. The complete trainer doesn't take up much room when stowed away. The adjustable holder for the simulated tissue stands firmly on the table without slipping, and the insert can be removed from the holder to be worked on while lying flat on the table.

Value-for-money tissue inserts can be purchased individually or in packs of 5.

Comes complete with: a silicone insert, a base, and a stand. Latex-free, silicone-based material.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorised distributor of 3B Scientific products in India with over 30 years of experience supplying medical education resources to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: What birth injuries can be practiced? A: Episiotomy (median, mediolateral, lateral), vaginal tears, and labial tears.

Q: Can the insert be used more than once? A: Yes — each cut can be sutured multiple times, and replacement inserts are available individually or in packs of 5.

Q: What's included? A: A silicone insert, a base, and a stand.

Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.