



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Life/form® Infant IV Arm-Injection & Punctures-1017949 [W44799]-LF03637U

Item No. LF03637U
Weight 0.7kg
Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories: Injections and Punctures



Product Description

A special, extremely thin, synthetic skin, and rubber tubing with appropriately small lumen and thin walls. The cephalic and basilic vein are accessible, as well as the dorsal venous arch on the hand. Comes with two IV bags with clamps, one pint Life/form® blood, and a winged infusion set. This product can be used on its own or added to the [Resusci Baby](#) (W19621 / 1017684) or [Life/form Crisis Baby](#) (W44717/1017248)

*Resusci® Baby is a trademark of the Laerdal Medical Corporation.

Why Choose the Life/form® Infant IV Arm for Your Nursing Skills Lab in India?

The Life/form® Infant IV Arm (Item No. LF03637U, Weight: 0.7kg, Make: Nasco Healthcare, USA) is built specifically for the scale and fragility of infant IV access — not a scaled-down adult arm. It uses an extremely thin synthetic skin along with rubber tubing sized with an appropriately small lumen and thin walls, so students learn to work with the delicacy real infant veins demand.

The cephalic and basilic veins are accessible, along with the dorsal venous arch on the hand — the key sites used for infant IV starts in clinical practice.

Flexible setup: This arm can be used on its own for focused IV-access practice, or mounted onto a Resusci® Baby or Life/form® Crisis Baby manikin for full-body infant resuscitation and emergency scenario training.

What's included: Two IV bags with clamps, one pint of Life/form® simulated blood, and a winged infusion set.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators and clinical skills training equipment to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: Which veins can be accessed on this arm? A: The cephalic and basilic veins, as well as the dorsal venous arch on the hand.

Q: Can this arm be used with a full manikin? A: Yes — it can be used standalone or added to a Resusci® Baby or Life/form® Crisis Baby manikin for full-scenario training.

Q: What's included with the arm? A: Two IV bags with clamps, one pint of Life/form® simulated blood, and a winged infusion set.

Q: How do I order from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.