



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



Geriatric LOR Kit (2) for Epidural and Spinal Injection Trainer

Item No. 1020629 [XP61-004]

Weight 1.66kg

Dimensions 28 x 14 x 75 cm

Brand 3B Scientific

[Read More](#)

SKU:

Categories: Injections and Punctures

Product Description

The geriatric LOR (Loss of Resistance) insert for the [Epidural and Spinal Injection Simulator P61](#) lets you practice with the specific anatomical deformations of the spine of an elderly patient. You can easily replace the standard patient LOR insert (delivered with your Epidural and Spinal Injection Simulator P61) with this optional geriatric insert. It represents the characteristic soft tissue structure (decreased resistance) and the deformed bone structure of a geriatric patient for advanced training possibilities:

Th11-Th12: No access (frustration threshold) Th12-L1: Centered, normal L1-L2: Ambilateral, small (partially bony center) L2-L3: Centered, small (dorsal process L3 shifted to left side) L3-L4: Centered, partially bony (dorsal process L3 shifted to left side) L4-L5: Centered, transverse oval Delivery content:

- 2 Geriatric Inserts
- 2 Spinal Tubes

Why Choose the Geriatric LOR Kit (2) for Your Epidural and Spinal Injection Simulator P61 in India?

The Geriatric LOR (Loss of Resistance) Kit — 2-pack (Item No. 1020629 [XP61-004], 1.66 kg, 28 × 14 × 7.5 cm*, Brand: 3B Scientific) gives institutions two geriatric inserts for the Epidural and Spinal Injection Simulator P61 — a cost-effective option for labs running higher training volumes or wanting a backup insert ready for continuous sessions.

Like the single-unit version, this insert replaces the standard patient LOR insert included with the P61, representing the characteristic soft tissue structure (decreased resistance) and deformed bone structure of a geriatric patient, including a deliberately non-accessible level (Th11-Th12) to build frustration-tolerance and technique recognition.

Vertebral-level detail:

- Th11-Th12: No access (frustration threshold)
- Th12-L1: Centered, normal
- L1-L2: Ambilateral, small (partially bony center)
- L2-L3: Centered, small (dorsal process L3 shifted to left side)
- L3-L4: Centered, partially bony (dorsal process L3 shifted to left side)
- L4-L5: Centered, transverse oval

Delivery contents: 2 geriatric inserts, 2 spinal tubes.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of 3B Scientific products in India with over 30 years of experience supplying medical simulators and clinical skills training accessories to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: How is this different from the single Geriatric LOR Insert? A: This kit includes two inserts and two spinal tubes instead of one — useful for higher training volumes or as a backup.

Q: Which simulator is this compatible with? A: The Epidural and Spinal Injection Simulator P61.

Q: How do I order from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.