



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Complete Infant CRiSis™ Manikin

Item No. LF03709U
Weight 8.204kg
Dimensions 61 x 42 x 19 cm
Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories:ALS Child

Product Description

Everything you need in one manikin! This dramatic, state-of-the-art training manikin is a complete resuscitation system designed for teaching life-saving techniques for infants. Ideal for pediatric advanced life support, nursing, hospital, paramedic, and EMT training at every level. Tactile and visual realism combine to provide students with the most realistic training possible. The most comprehensive PALS trainer available, the manikin allows practice of 4-lead ECG monitoring, IV therapy, umbilical catheterization, airway management, intraosseous infusion, and CPR - all in one manikin! The unique chest skin with an umbilical cord allows practice of proper techniques in umbilical catheterization, including cutting and clamping. Fluid may be aspirated through the umbilicus, making this training manikin incredibly lifelike. An IV arm and leg allow for realistic practice of pediatric venipuncture. The cephalic and basilic veins are accessible on the IV arm, as well as the dorsal venous arch on the hand; the great saphenous veins are accessible on the IV leg, as well as the dorsal venous arch on the foot.

An intraosseous infusion leg demonstrates and simulates the intraosseous infusion procedure. Featuring realistic anatomy and size, the airway management head is ideal for practicing airway management skills on an infant. Includes a hard carrying/storage case. Procedures That Can Be Performed: CPR: Palpable and visual landmarks; fully articulated head, neck, and jaw; jaw thrust Airway Management: Realistic anatomy of the mouth, tongue, oral pharynx, larynx, epiglottis, vocal cords, trachea, and esophagus; oral and digital intubation capabilities; suctioning capabilities IV Arm & Leg: Smaller lumen tubing simulates the size of an infant's veins; realistic flashback from a pressurized system; replaceable skin and veins ECG/Umbilical Cannulation Skin Arrhythmia recognition; learn proper techniques in umbilical catheterization; allows practice of umbilical cannulation, cutting, and clamping; aspirate fluid through the umbilicus Intraosseous Infusion: Palpable landmarks include the patella, tibia, and tibial tuberosity; replaceable bones and skin; pressurized system allows aspiration of fluid

Why Choose the Life/form® Complete Infant CRiSis™ Manikin for Your Simulation Lab in India?

The Life/form® Complete Infant CRiSis™ Manikin (Item No. LF03709U / 1005658) by Nasco Healthcare, USA is the most comprehensively equipped configuration in the Infant CRiSis™ family that does not include the Interactive ECG Simulator as standard — 8.204 kg, 61 × 42 × 19 cm, with a hard carry case included. It adds the ECG/Umbilical Cannulation Skin over the Intermediate Infant CRiSis™ (1021997), making it the first configuration in the Infant family where umbilical catheterisation training — a critical neonatal vascular access skill unique to this age group — is available alongside airway management, CPR, IO infusion, IV arm, and IV leg. For NMC-mandated simulation labs, neonatal resuscitation training centres, NICU in-service training departments, and paediatric nursing programmes in India that need a single infant manikin covering the complete NRP and PALS neonatal procedural skill set, this is the correct configuration.

The ECG/Umbilical Cannulation Skin is this manikin's defining capability. The chest skin features a realistic umbilical cord stump that enables trainees to practise the full umbilical catheterisation sequence — cutting, clamping, and cannulation — with fluid aspiration through the umbilicus confirming correct technique. This skill cannot be practised on any manikin without this specialised skin, and the Complete configuration is the entry-level point where it becomes available. The IV arm provides access to cephalic and basilic veins and the dorsal venous arch on the hand; the IV leg provides access to the great saphenous veins and the dorsal venous arch on the foot — both using smaller lumen tubing that simulates the size of an infant's veins with realistic pressurised flashback and replaceable skin and veins. The IO infusion leg features palpable patella, tibia, and tibial tuberosity landmarks with a pressurised aspiration system and replaceable bones and skin. The airway management head presents realistic infant anatomy — mouth, tongue, oral pharynx, larynx, epiglottis, vocal cords, trachea, and oesophagus — and supports oral and digital intubation and suctioning. CPR features palpable and visual landmarks, fully articulated head, neck, and jaw, and jaw thrust. The Interactive ECG Arrhythmia Simulator (1005666) is optional for this configuration and included as standard only in the Deluxe Infant CRiSis™ (1018146).

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators to 1,500+ institutions across 28 states. We provide pan-India delivery, on-site installation, faculty orientation, and dedicated after-sales support. All products carry genuine manufacturer warranty and meet ISO 9001 quality standards.

Frequently Asked Questions

Q: What does this manikin include that the Intermediate Infant CRiSis™ (1021997) does not? A: The ECG/Umbilical Cannulation Skin — enabling umbilical catheterisation training including cutting, clamping, cannulation, and fluid aspiration through the umbilicus. The Intermediate has IO infusion, IV arm, and IV leg but does not have this skin.

Q: Is the Interactive ECG Arrhythmia Simulator included? A: No. The Interactive ECG Simulator (1005666) is an optional add-on for this configuration. It is included as standard only in the Deluxe Infant CRiSis™ (1018146).

Q: What is the warranty and how do I get support? A: All Nasco Healthcare products carry standard manufacturer warranty. For support, reach us at sem@semtrainers.com or call +91-88495 63724.

Manikin can be used with Supraglottic airway devices for airway training.

| | Basic Infant CRiSis™ | Intermediate Infant CRiSis™ | Complete Infant CRiSis™ | Deluxe Infant CRiSis™ |
|-------------------|-------------------------|-----------------------------|-------------------------|-------------------------|
| Item No. | 1017248 | 1021997 | 1005658 | 1018146 |
| Airway management | ✓ | ✓ | ✓ | ✓ |

| | | | | |
|---|--------------------------------------|--------------------------------------|--------------------------------------|----------|
| CPR | ✓ | ✓ | ✓ | ✓ |
| ECG/Umbilical cannulation Skin | - | - | ✓ | ✓ |
| Intraosseous infusion | - | ✓ | ✓ | ✓ |
| IV arm | - | ✓ | ✓ | ✓ |
| IV leg | - | ✓ | ✓ | ✓ |
| Interactive ECG arrhythmia simulator | optional (1005666) | optional (1005666) | optional (1005666) | included |
| Carry case | included | included | included | included |