



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Smart Stat Xtreme Trauma Bleeding Leg

Item No. 101-482
Weight 10 lbs.
Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories: Moulage and Wound Simulation

Product Description

These graphic explosion wounds offer detailed enhancements to challenge any practitioner's level of care. These moulages feature missing muscle and bone, flesh avulsions, fatty tissue, and bone ends. The blood supply line from the SMART STAT simulator connects quickly, without blood splatter, to the moulage, and has pressurized blood flow to create a messy wound. Proper application of a tourniquet will halt the blood flow. Xtreme Bleeding Moulages can be used with the STAT line of patient simulators. External blood supplies must be used if not attached to the SMART STAT. Consists of a molded moulage piece and an interior leg structure. This feature provides a simple connection to the existing knee joint of the SMART STAT patient simulator using an easy-to-install knee pin. Ship weight 10 lbs.

Smart Stat Xtreme Trauma Bleeding Leg (Nasco Healthcare) — Pressurized

Bleeding Moulage for STAT Simulators in India

The Smart Stat Xtreme Trauma Bleeding Leg (Item No. 101-482, Brand: Nasco Healthcare, USA) is a graphic explosion-wound moulage designed to challenge advanced levels of trauma care, featuring missing muscle and bone, flesh avulsions, fatty tissue, and visible bone ends.

The blood supply line from the SMART STAT simulator connects quickly to the moulage without blood splatter, delivering pressurized blood flow to create a realistic, messy wound. A properly applied tourniquet will halt the blood flow, giving trainees a genuine feedback loop for hemorrhage control technique. The moulage consists of a molded piece and an interior leg structure that connects to the existing knee joint of the SMART STAT patient simulator via an easy-to-install knee pin. An external blood supply is required if the moulage is not attached directly to a SMART STAT unit.

Specifications: Weight (ship) 10 lbs (approx. 4.54 kg). Brand: Nasco Healthcare, USA. Compatible with the STAT line of patient simulators.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying advanced trauma and hemorrhage-control training equipment to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

FAQ Q: How does the pressurized bleeding feature work? **A:** The SMART STAT blood supply line connects to the moulage without splatter, delivering pressurized flow that a correctly applied tourniquet will stop — giving trainees real feedback on hemorrhage control technique.

Q: What simulators is this compatible with? **A:** The STAT line of Nasco patient simulators, connecting via an easy-to-install knee pin at the existing knee joint.

Q: How do I order from SEM Trainers? **A:** Contact sem@semtrainers.com or call +91-88495 63724 for pricing and availability.