



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com

Life/form® Advanced Injection Belly

Item No. LF09900


Weight

Make: Nasco Healthcare, USA

[Read More](#)

SKU:LF09900

Categories:Injections and Punctures



Product Description

The Life/form® Advanced Injection Belly is an invaluable tool for showing how to properly insert and rotate the infusion set for insulin pumps, avoiding the 2 in. (5 cm) area around the navel, as well as giving many types of self-injections (including pen-style injections). Similar to the Life/form® Injection Belly, but made of a material that creates a more realistic experience. The trainer has lifelike skin with what feels like real stomach tissue, along with the ability to pinch the skin. Backing keeps needles from going through. Adjustable straps allow the stomach to be attached to a manikin or human patient for simulation practice. The soft material can be punctured repeatedly without showing holes. Liquid injection is not recommended. Dimension: 10 in. L x 6 in. W x 2 in. H. Manufactured using a quality management system certified as being in conformity with ISO 9001.

Good — solid content here. The genuinely distinctive angle: this is upgraded material over the standard

Injection Belly, specifically built for **insulin pump infusion set rotation practice** (a skill most injection trainers don't address) plus general self-injection/pen-style injection practice — and it explicitly calls out that **liquid injection is not recommended**, which buyers expecting a fluid-capable trainer should know upfront.

Why Choose the Life/form® Advanced Injection Belly for Your Nursing Skills Lab in India?

The Life/form® Advanced Injection Belly (Item No. LF09900, Make: Nasco Healthcare, USA) is built for teaching proper insulin pump infusion set insertion and site rotation — including how to avoid the 2 in. (5 cm) zone around the navel — as well as general self-injection technique, including pen-style injections.

It's an upgraded version of the standard Life/form® Injection Belly, made from a material that delivers a more realistic feel: lifelike skin over what feels like real stomach tissue, with the ability to pinch the skin as you would on a real patient. A firm backing keeps needles from passing through, and the soft material can be punctured repeatedly without showing visible holes.

Note: Liquid injection is not recommended with this trainer — it's designed for needle insertion and site-rotation technique, not fluid delivery practice.

Adjustable straps let the belly attach to a manikin or a human patient for hands-on simulation practice.

Specifications:

- Dimensions: 10" L × 6" W × 2" H
- Manufactured under a quality management system certified to ISO 9001

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators and clinical skills training equipment to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: Can I use liquid injections with this trainer? A: No — liquid injection is not recommended. It's designed for needle insertion and infusion-set rotation practice.

Q: Can this belly be attached to a manikin or a human patient? A: Yes — adjustable straps allow it to be worn on a manikin or a human patient for realistic simulation practice.

Q: How is this different from the standard Life/form® Injection Belly? A: It uses a more advanced

material that gives a more realistic tissue feel, including the ability to pinch the skin.

Q: How do I order from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.