

**SEM****SEM TRAINERS & SYSTEMS**

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Life/form® Advanced Suture Kit

Item No. LF00894
Weight
Make: Nasco Healthcare, USA

[Read More](#)

SKU:LF00894

Categories:Suturing and Bandaging

Product Description

The Advanced Suture Kit offers a soft, flexible three-layer skin opens realistically when incised, yet is highly resistant to tearing - sutures won't cut through even when pulled tight. Skin pad can be sutured hundreds of times. Includes easy-to-change specialty inserts for deep dermal suturing, tendon repair, and purse-string suture - challenging procedures that require a high degree of precision. Compact, portable and completely washable. Includes base, one skin pad, four non-latex bands, one muscle block, one deep dermal insert, one tendon block, one tendon, one purse-string insert, suture kit, and carry bag.

Item No. LF00894, weight field blank, Make: Nasco Healthcare, USA ?. This is the most advanced suture trainer in the subcategory — three specialty inserts for procedures that simple skin pad suturing can't train: deep dermal suturing, tendon repair, and purse-string suture. Three-layer skin pad that

opens realistically when incised but can be sutured hundreds of times without tearing. Downloadable manual available.

DESCRIPTION (copy-paste into WordPress product description):

Why Choose the Life/form® Advanced Suture Kit for Your Simulation Lab in India?

The Life/form® Advanced Suture Kit (Item No. LF00894) by Nasco Healthcare, USA is a comprehensive, portable suturing skills trainer that goes significantly beyond standard skin pad suturing — three interchangeable specialty inserts train deep dermal suturing, tendon repair, and purse-string suture technique, all procedures requiring a degree of precision that standard flat skin pad trainers cannot provide. For surgical skills training programmes, medical college simulation labs, emergency medicine training, and NMC-mandated clinical skills labs in India that need to train the full range of surgical suturing techniques beyond basic interrupted and continuous skin closure, this kit provides the complete advanced suturing curriculum on a single compact, washable platform. A downloadable user manual is available on the product page.

The three-layer skin pad is the foundation of the kit — soft, flexible, and realistic in feel, it opens cleanly when incised and yet is highly resistant to tearing, so sutures won't cut through even when pulled tight during practice. The skin pad can be sutured hundreds of times, making the cost per training session exceptionally low compared to single-use pad systems. Four non-latex bands hold the pad securely on the base during practice. The three specialty inserts convert the same base platform into three distinct training stations: the deep dermal insert trains suturing below the skin surface where suture placement cannot be directly visualised; the tendon block with included tendon trains the precision alignment and tensioning required for tendon repair technique; and the purse-string insert trains the circular closure technique used in bowel anastomosis, appendiceal stump closure, and port site closure. Compact and completely washable, the kit is suitable for both lab and outreach training.

Includes: base, one skin pad, four non-latex bands, one muscle block, one deep dermal insert, one tendon block, one tendon, one purse-string insert, suture kit, and carry bag.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators to 1,500+ institutions across 28 states. We provide pan-India delivery, on-site installation, faculty orientation, and dedicated after-sales support.

Q: What advanced suturing techniques does this kit train? A: Deep dermal suturing, tendon repair, and purse-string suture — three procedures requiring precision that cannot be trained on a standard flat skin pad, in addition to standard skin closure on the three-layer skin pad.

Q: How many times can the skin pad be sutured? A: The skin pad can be sutured hundreds of times and is highly resistant to tearing — sutures will not cut through even when pulled tight.

Q: What is included? A: Base, one skin pad, four non-latex bands, one muscle block, one deep dermal insert, one tendon block, one tendon, one purse-string insert, a suture kit, and a carry bag.