



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



Nose Bleed Trainer

Item No. SB38892

Weight

Make: Nasco Healthcare, USA

[Read More](#)

SKU:SB38892

Categories:Nose and Throat Examination

Product Description

- Interactive model provides a practical test of skill in managing severe nose bleeds
- Students must identify the bleeding point, clear it by suction, and stop the bleeding either by nasal packing or cautery
- The nasal interior has capillary tubes for delivering artificial blood at variable rates to four sites on the septum
- This Interactive model also contains an indicator light that illuminates to signal a successful cautery
- Includes hard plastic head with soft nasal area, control box, four 10 ml syringes, a 50 ml syringe, blood concentrate, and instruction manual
- Requires one 9V battery (included)

- 11 in. L x 6-1/3 in. W x 7-1/2 in. H
- Two-year warranty

Excellent existing content — one of the best bullet-point descriptions in this subcategory. Item No. SB38892, weight field blank, Make: Nasco Healthcare, USA ?. Category correctly placed under Nose and Throat Examination only. Size 11" x 6 2/3" x 7 1/2" (approx. 28 x 16 x 19 cm). Four capillary tube bleeding sites on the septum, variable flow rates, cautery success indicator light, nasal packing training. Two-year warranty. 9V battery included.

DESCRIPTION (copy-paste into WordPress product description):

Why Choose the Nose Bleed Trainer for Your Simulation Lab in India?

The Nose Bleed Trainer (Item No. SB38892) by Nasco Healthcare, USA is an interactive epistaxis management trainer — a hard plastic head with soft nasal area (approx. 28 x 16 x 19 cm) that provides a practical test of skill in managing severe nose bleeds by requiring trainees to identify the bleeding point, clear it by suction, and stop the bleeding by nasal packing or cautery. For ENT residency programmes, emergency medicine training, general practice training, and NMC-mandated clinical skills labs in India that need to train active epistaxis management rather than just anatomy, this trainer provides the full clinical skill sequence — identification, suction clearance, and haemostasis — on a single interactive model. Backed by a two-year warranty.

The nasal interior has four capillary tubes delivering artificial blood at variable rates to four distinct sites on the septum, allowing the instructor to simulate different epistaxis presentations and bleeding severities during a single session. Trainees must use suction to clear the simulated blood to identify which of the four septal sites is bleeding — training the active clearance step that precedes both packing and cautery in clinical epistaxis management. The model includes an indicator light that illuminates to signal a successful cautery, providing objective confirmation that the trainee has correctly performed the cautery technique on the correct site. Nasal packing technique can also be practised on the soft nasal area.

Includes: hard plastic head with soft nasal area, control box, four 10 ml syringes, one 50 ml syringe, blood concentrate, instruction manual, and one 9V battery.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators to 1,500+ institutions across 28 states. We provide pan-India delivery, on-site installation, faculty orientation, and dedicated after-sales support.

Frequently Asked Questions

Q: How many bleeding sites does this trainer have? A: Four capillary tube sites on the septum, each delivering artificial blood at variable rates — allowing the instructor to simulate different epistaxis presentations during a session.

Q: Does the trainer provide feedback on successful cautery? A: Yes. An indicator light illuminates when cautery has been correctly performed on the active bleeding site.

Q: What is included and what is the warranty? A: Hard plastic head with soft nasal area, control box, four 10 ml syringes, one 50 ml syringe, blood concentrate, instruction manual, and a 9V battery. Two-year warranty.