

**SEM****SEM TRAINERS & SYSTEMS**

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724**Mobile 1:** +91 98791 03905**Email:** sem@semtrainers.com**Website:** www.semtrainers.com

Suprapubic Training Model - Light

Item No. SB38897 L**Weight****Make:** Nasco Healthcare, USA[Read More](#)**SKU:**SB38897**Categories:**Catheterization

Product Description

Become skilled in this increasingly common alternative to urethral indwelling catheters. Feel the normal pressures exerted by the abdominal wall on catheterization. Understand the average length of catheter required to be inserted in order to drain the bladder. Water will flow through catheter when bladder is successfully reached. Lower torso is attached to a vacuum-formed base. Also includes plastic container with lid and fittings, pressure bulb with tubing and connectors, size 16 Foley catheter, instruction manual, and rigid carrying case.

A genuine Nasco counterpart to the Kyoto Kagaku suprapubic products covered just before this — same skill, different brand, worth positioning as an alternative option. One data gap (missing weight) and one naming pattern worth flagging based on precedent from earlier in this session.

[Cert] Item No. SB38897 L, Make Nasco Healthcare, USA. No weight given. [Guessing] The "L" suffix strongly suggests this is a Light skin-tone variant, following the same naming convention as the Male/Female Catheterization Simulator (LF00843 Light / LF00842 Medium) covered earlier this session — a Medium variant (likely SB38897 M or similar) may exist separately. Worth confirming with your supplier before publishing, since a buyer wanting Medium skin tone would need a different item number entirely.

[Cert] This trains suprapubic catheterization specifically — described as "an increasingly common alternative to urethral indwelling catheters" — with realistic abdominal wall pressure feedback and a fluid-flow confirmation when the bladder is successfully reached, teaching students the average catheter insertion length required.

DESCRIPTION (copy-paste into WordPress):

Why Choose the Suprapubic Training Model — Light — for Your Nursing Skills Lab in India?

The Suprapubic Training Model, Light (Item No. SB38897 L, Make: Nasco Healthcare, USA) builds skill in suprapubic catheterization — an increasingly common alternative to urethral indwelling catheters, particularly relevant for patients who cannot be catheterized transurethrally.

Students feel the normal pressures the abdominal wall exerts during catheterization and learn the average catheter insertion length required to successfully drain the bladder. Fluid flows through the catheter once the bladder is correctly reached, giving clear, unambiguous confirmation of successful placement. The lower torso is mounted on a vacuum-formed base for stable table-top training.

Note on skin tone: This listing is the Light skin tone version (SB38897 L). Please confirm with SEM Trainers whether a Medium skin tone variant is available under a different item number before ordering.

Includes: Plastic container with lid and fittings, pressure bulb with tubing and connectors, size 16 Foley catheter, instruction manual, rigid carrying case.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators and clinical skills training equipment to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: What catheter size does this use? A: A size 16 Foley catheter (included).

Q: Is a Medium skin tone available? A: Please contact SEM Trainers to confirm availability of alternative skin tone options.

Q: How do I order from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.