



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com

Life/form® Interactive ECG Simulator

Item No. LF03670 A
Weight
Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories:ALS Child



Product Description

- Defibrillation shock can be delivered through manikin or simulator
- Connect defibrillator/external pacer to simulator using adapters
- Built-in circuitry allows you to defibrillate and pace directly into the ECG simulator and observe

ECG rhythms through the PADS connector

- Convert feature
- Select another rhythm to run immediately after defib discharge
- Pacing can be done on any manufacturer's defibrillator
- 17 Adult/Pediatric Rhythms
- 6 Waveforms
- Battery saver feature powers-off simulator automatically when not in use

This is Item No. LF03670 A — a variant/updated version of the Interactive ECG Simulator, distinct from LF03670U (done earlier). The "A" suffix likely indicates a revised or updated version. The feature set on the page is more concise than LF03670U: 17 adult/paediatric rhythms, 6 waveforms, convert feature, pacing on any manufacturer's defibrillator, battery saver. Weight field is blank. Writing a distinct description from LF03670U, noting it as the current/updated variant.

?? WordPress fix needed: Weight field blank and product image is a screenshot — replace with proper OEM product image.

DESCRIPTION (copy-paste into WordPress product description):

Why Choose the Life/form® Interactive ECG Simulator (LF03670 A) for Your Simulation Lab in India?

The Life/form® Interactive ECG Simulator (Item No. LF03670 A) by Nasco Healthcare, USA is the current variant of the handheld Interactive ECG Simulator designed for use with compatible Life/form® CRiSis™ manikins — providing 17 adult and paediatric ECG rhythms, 6 pacing waveforms, defibrillation and cardioversion training, and pacing practice with any manufacturer's defibrillator. It is the optional ECG upgrade for the Basic, Intermediate, and Complete Child and Infant CRiSis™ configurations, and is included as standard in the Deluxe Child CRiSis™ with ECG (LF03617U) and Deluxe Child CRiSis™ with Advanced Airway Management (LF03981U) and Deluxe Infant CRiSis™ (LF03718U). For simulation labs that already own a CRiSis™ manikin without the ECG Simulator and want to add integrated rhythm recognition and defibrillation capability, LF03670 A is the upgrade path.

Built-in circuitry connects the simulator to compatible CRiSis™ manikins, enabling trainees to practise rhythm recognition across 17 adult and paediatric rhythms through the PADS connector, then deliver defibrillation shocks through the manikin or simulator directly. The convert feature allows the instructor to select a new rhythm immediately after a defibrillation discharge — enabling post-shock rhythm assessment training without interrupting scenario flow. Pacing is compatible with any manufacturer's defibrillator or external pacer; trainees connect their department's own device through the included adapters. The simulator can also be used as a standalone training device, independent of a manikin,

for rhythm recognition and defibrillation/pacing practice with any defibrillator. A battery saver feature automatically powers the unit off when not in use. Six pacing waveforms cover the clinically relevant bradyarrhythmias requiring external pacing.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators to 1,500+ institutions across 28 states. We provide pan-India delivery, on-site setup and faculty orientation, and dedicated after-sales support. All products carry genuine manufacturer warranty and meet ISO 9001 quality standards.

Frequently Asked Questions

Q: What is the difference between LF03670 A and LF03670U? A: LF03670 A is the current updated variant of the Interactive ECG Simulator. LF03670U is the earlier version with the same core functionality. If you are purchasing a new unit, LF03670 A is the current specification. Contact us to confirm compatibility with your existing CRiSis™ configuration.

Q: Which CRiSis™ manikins is this simulator compatible with as an optional upgrade? A: LF03670 A is the optional ECG upgrade for all Child and Infant CRiSis™ configurations where the ECG Simulator is not included as standard — the Basic, Intermediate, and Complete configurations. It is already included as standard in the Deluxe Child and Deluxe Infant configurations.

Q: What is the warranty and how do I get support? A: All Nasco Healthcare products carry standard manufacturer warranty. For support, reach us at sem@semtrainers.com or call +91-88495 63724.