



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com





Daily QA Phantom WEM "Aqua Slab"

Item No. PH-80
Weight
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories: Radiology Phantom



Product Description

Features

Training skills / Applications

Phantom using innovative Water Equivalent Material Aqua Slab Save time and reduce risk with relate to water phantoms
Daily QA of CT

Size (approx.)
Weight (approx.)
Materials
Related products
Update

41948-000:W20 x D20 x H20 cm	41948-100:W20 x D20 x H21.5 cm
41948-000 / 6.4kg	41948-100 / 6.5kg
Polyurethane (Aqua Slab)	
PH-75A Multi Energy CT Quality Assurance Phantom (TR-J) PH-75B Multi Energy CT Quality Assurance Phantom (TR-J)	
September 1, 2020	

Why Choose the Daily QA Phantom WEM "Aqua Slab" for CT Quality Assurance in India?

The Daily QA Phantom WEM "Aqua Slab" (Item No. PH-80, Brand: Kyoto Kagaku, Japan) uses innovative Water Equivalent Material (Aqua Slab) construction to save time and reduce the risks associated with traditional water phantoms — a practical daily QA solution for busy CT departments.

Available in two sizes to match different CT QA workflow needs.

Training applications:

- Daily QA of CT systems

Specifications:

- 41948-000: W20 × D20 × H20cm — Weight approx. 6.4kg
- 41948-100: W20 × D20 × H21.5cm — Weight approx. 6.5kg
- Materials: Polyurethane (Aqua Slab)

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Kyoto Kagaku products in India with over 30 years of experience supplying radiology QA phantoms and CT quality assurance equipment to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: What is Aqua Slab material, and why use it instead of a water phantom? A: Aqua Slab is a Water Equivalent Material in solid polyurethane form — it saves setup time and reduces the handling risks associated with traditional liquid-filled water phantoms while maintaining equivalent QA performance.

Q: Which size should I choose? A: Both sizes (41948-000 and 41948-100) support daily CT QA; the choice depends on your specific QA protocol requirements — contact us for guidance based on your CT system.

Q: How do I order from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.