



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE


SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



Acrylic Phantom XAC

Item No. PH-14

Weight

Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories: Radiology Phantom

Product Description

Slab phantoms for radiation absorption and scattering measurement

Training skills / Applications

Set includes

Materials

X-ray

1 phantom / storage case

Acrylic resin

41430-000	XAC-01	30x30x0.1cm/11.8x11.8x0.04inch	41431-000
XAC-02	30x30x0.2cm/11.8x11.8x0.08inch	41432-000	XAC-03
30x30x0.3cm/11.8x11.8x0.12inch	41433-000	XAC-04	
30x30x0.4cm/11.8x11.8x0.16inch	41434-000	XAC-05	
30x30x0.5cm/11.8x11.8x0.2inch	41435-000	XAC-08	
30x30x0.8cm/11.8x11.8x0.3inch	41436-000	XAC-1	
30x30x1cm/11.8x11.8x0.4inch	41437-000	XAC-2	
30x30x2cm/11.8x11.8x0.8inch	41438-000	XAC-3	
30x30x3cm/11.8x11.8x1.2inch	41439-000	XAC-4	
30x30x4cm/11.8x11.8x1.6inch	41440-000	XAC-5	
30x30x5cm/11.8x11.8x2.0inch	41441-000	XAC-8	
30x30x8cm/11.8x11.8x3.1inch	41442-000	XAC-10	
30x30x10cm/11.8x11.8x3.9inch	*with storage case		
July 10, 2020			

other versions

Update

Acrylic Phantom XAC (PH-14) — Slab Phantoms for X-ray Absorption & Scattering Measurement

The Acrylic Phantom XAC (Item No. PH-14) by Kyoto Kagaku, Japan is a series of acrylic slab phantoms built for measuring radiation absorption and scattering in X-ray applications — a standard tool for radiography quality assurance and dosimetry testing.

13 thickness options for precise, incremental testing. The XAC series spans from 0.1cm to 10cm thick, all at a consistent 30 × 30 cm footprint: XAC-01 (0.1cm), XAC-02 (0.2cm), XAC-03 (0.3cm), XAC-04 (0.4cm), XAC-05 (0.5cm), XAC-08 (0.8cm), XAC-1 (1cm), XAC-2 (2cm), XAC-3 (3cm), XAC-4 (4cm), XAC-5 (5cm), XAC-8 (8cm), and XAC-10 (10cm). This range lets institutions build custom test stacks by combining slabs to reach specific total thicknesses, or run graduated absorption/scattering measurements across a controlled range.

Simple, purpose-built construction. Made from acrylic resin, each phantom is supplied with a storage case for safe handling and storage between uses.

Why this matters for radiography QA programs: consistent, known-thickness reference phantoms are essential for calibrating and verifying X-ray equipment performance, particularly for absorption and scatter characteristics that affect image quality and radiation dose.

What's included (per unit): 1 phantom, storage case.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorised distributor of Kyoto Kagaku products in India with over 30 years of experience supplying medical simulators and clinical skills trainers to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: Which thickness should I order? A: It depends on your specific QA or research testing needs — the series spans 0.1cm to 10cm. Contact our team with your target measurement range and we'll help identify the right thickness (or combination of thicknesses) for your application.

Q: Can multiple slabs be combined? A: Yes — since all slabs share the same 30 × 30 cm footprint, they can be stacked to reach custom total thicknesses for specific test scenarios.

Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability, specifying which thickness(es) you need.