



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Water Body Phantom WAC

Item No. PH-17
Weight 9 kg
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories: Radiology Phantom

Product Description

Image evaluation in plain X-ray

Features

Training skills / Applications

Set includes

Size (approx.)

Packing size (approx.)

Packing weight (approx.)

Materials

Update

Water body phantom represents human chest and abdomen to serve as radiation absorber and scatterer.
X-ray
1 body phantom / 1 storage case
W30 x D20 x H45 cm W11.8 x D7.9 x H17.7 in
W57x D34xH44 cm / W22.4xD13.4xH17.3 inch
9kg / 19.8lb
Acrylic resin
July 10, 2020

Why Choose the Water Body Phantom WAC for Plain X-Ray Image Evaluation in India?

The Water Body Phantom WAC (Item No. PH-17, Brand: Kyoto Kagaku, Japan) represents the human chest and abdomen, serving as a radiation absorber and scatterer for plain X-ray image evaluation training.

Training applications:

- X-ray image evaluation

Set includes: 1 body phantom, 1 storage case.

Specifications:

- Size approx. W30×D20×H45cm (W11.8×D7.9×H17.7in)
- Packing size approx. W57×D34×H44cm (W22.4×D13.4×H17.3in)
- Packing weight approx. 9kg (19.8lb)
- Materials: Acrylic resin

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Kyoto Kagaku products in India with over 30 years of experience supplying plain X-ray training and QA equipment to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: What is this phantom used for? A: It serves as a radiation absorber and scatterer, simulating human chest and abdomen tissue for evaluating plain X-ray image quality.

Q: How do I order from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.