



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com




Rectal and Bladder Ultrasound Phantom

Item No. US-23
Weight
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories:Ultrasound Phantom



Product Description

Ultrasound training cubes for assessment of residual urine and excretion

Features

Training skills / Applications

Set includes

Size (approx.)

Weight (approx.)

Update

Variation of residue excretion in the rectum for ultrasound assessment: hard, soft, combination and empty
Three types of urinal bladder for assessment and urine volume measurement: small, middle, large with a balloon catheter

Ultrasound assessment of residue excretion in the rectum
Measurement of urine volume
Ultrasound assessment of the urinal bladder and the prostate
Assessment of indwelling bladder catheter

3 types of ultrasound phantom cube / 1 talcum powder / 1 storage case / 1 instruction manual

W15 x D15 x H20 cm

2.0 kg (each)

March 14, 2024

Rectal and Bladder Ultrasound Phantom (US-23) — Residual Urine & Bowel Content Assessment Training

The Rectal and Bladder Ultrasound Phantom (Item No. US-23) by Kyoto Kagaku, Japan is a set of ultrasound training cubes built specifically for assessing residual urine volume and rectal excretion content by ultrasound — two closely related bedside assessment skills often taught together in nursing and geriatric care training.

Rectal assessment variations. The rectal training cube includes variations in residue consistency — hard, soft, combination, and empty — letting students practice distinguishing different presentations by ultrasound rather than training on a single fixed appearance.

Bladder assessment across three sizes. Three urinary bladder variants (small, middle, large), each with a balloon catheter, support both urine volume measurement and assessment of an indwelling bladder catheter — a common bedside nursing skill for monitoring patients with catheters in place.

Training skills covered: ultrasound assessment of rectal residue, urine volume measurement, ultrasound assessment of the urinary bladder and prostate, and assessment of indwelling bladder catheters.

What's included: 3 types of ultrasound phantom cube, talcum powder, storage case, instruction manual. **Size:** W15 × D15 × H20 cm. **Weight:** 2.0 kg (each cube).

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorised distributor of Kyoto Kagaku products in India with over 30 years of experience supplying medical simulators and clinical skills trainers to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: Can this be used to train catheter care, or only bladder scanning? A: Both — the balloon catheter on each bladder phantom supports assessment of an indwelling catheter, alongside urine volume measurement.

Q: What does the talcum powder in the kit get used for? A: It's typically used as a lubricant/gel substitute or surface prep aid for ultrasound probe contact during training — contact our team if you'd like the manufacturer's specific usage guidance.

Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.