



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com





Scrotal Ultrasound Phantom

Item No. US-11
Weight 7 kg
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories:Ultrasound Phantom



Product Description

Master techniques for a safe and painless examination

Features

1. Two types of interchangeable scrotums: normal and pathology
2. Excellent ultrasound image quality

Training skills / Applications

Case / Pathology

Set includes

Size (approx.)

Packing size (approx.)

Weight (approx.)

Packing weight (approx.)

Materials

Replacement parts

Related products

Production ? Development Supervision

Update

Patient positioning and preparation for examination / Scrotal ultrasound screening / Visualization of testicular cancer
Normal Pathological
1 pelvic body / 1 normal scrotum / 1 pathological scrotum / 1 talcum powder / 1 storage case / 1 instruction manual
W34 x D33 x H24cm / W13.4 x D13 x H9.4 in
W52 x D45 x H38cm / W20 x D18 x H15in
3.6 kg / 7.9 lbs
7 kg / 15.4 lb
Soft resin
41911-010 normal scrotum unit 41911-020 pathological scrotum unit
US-10A Female Pelvic Ultrasound Phantom Full set US-16 Intravesical Urine Volume Measurement Phantom
Charlotte Henningsen, MS, RT(R), RDMS, RVT, FSDMS, FAIUM Chair & Professor - Sonography Department Adventist University of Health Sciences
August 21, 2020

Scrotal Ultrasound Phantom (US-11) — Screening & Testicular Pathology Training

The Scrotal Ultrasound Phantom (Item No. US-11, 3.6 kg / packing weight 7 kg) by Kyoto Kagaku, Japan is built to help trainees master safe, patient-comfortable technique for scrotal ultrasound examination — a procedure where positioning and handling matter as much as image interpretation.

Two interchangeable scrotums, normal and pathological, let students practice on both a healthy baseline and a pathological presentation without needing two separate phantom units. This directly supports one of the most clinically important applications of scrotal ultrasound: visualization of testicular cancer, alongside general screening and patient positioning/preparation technique.

Excellent ultrasound image quality ensures the training experience closely mirrors real clinical scanning, so skills transfer effectively to actual patient examinations.

What's included: 1 pelvic body, 1 normal scrotum, 1 pathological scrotum, talcum powder, storage case, instruction manual. **Size:** W34 × D33 × H24 cm. Made from soft resin.

Replacement parts available for both the normal scrotum unit and pathological scrotum unit separately — useful if one wears out faster than the other under heavy training use, letting you replace just what's needed rather than the full kit.

Developed under supervision of Charlotte Henningsen, Chair & Professor of the Sonography

Department at Adventist University of Health Sciences.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorised distributor of Kyoto Kagaku products in India with over 30 years of experience supplying medical simulators and clinical skills trainers to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: Can students practice on both normal and pathological anatomy with one unit? A: Yes — the kit includes both a normal and a pathological scrotum, interchangeable on the same pelvic body.

Q: What pathology does the pathological scrotum represent? A: It's built to support visualization training for testicular cancer, alongside general pathological screening technique.

Q: Are individual scrotum units replaceable? A: Yes — the normal and pathological scrotum units are each available as separate replacement parts.

Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.