



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com


Infant Hip Sonography Training Phantom

Item No. US-13
Weight 7.7 kg
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories:Ultrasound Phantom



Product Description

World's first infant hip sonography phantom/Best tool to teach Graf's method

Features

1. The market's only training phantom for hip sonography on a full body manikin of 6-week-old infant
2. Bilateral hips for examination
3. Facilitates anatomical understanding

Training skills / Applications

Set includes

Size (approx.)

Packing size (approx.)

Weight (approx.)

Packing weight (approx.)

Materials

Related products

Production ? Development Supervision

Update

Setting and preparation for hip sonography / Changing the position of the infant / Communication and interaction with infant's guardian / Correct use and positioning of the transducer / Recognition of landmarks for hip sonography / Visualization of standard, anterior and posterior planes / Interpretation and morphological classification of the sonogram

1 ultrasound infant phantom / 1 talcum powder / 1 storage case / 1 instruction manual

W55 x D25 x H13 cm W21.6 x D9.8 x H5.1 in

W71 x D36 x H22cm / W28.0 x D14.2 x H8.7in

3kg / 6.6lb

7.7kg / 17lbs

Polyurethane elastomer/ Latex free

[US-14a Ultrasound Neonatal Head Phantom \(Normal type\)](#)

[US-14b Ultrasound Neonatal Head Phantom \(Abnormal type\)](#)

Univ. Prof., h.c., Reinhard Graf, M.D. Japanese Society of Orthopedic Ultrasonics

August 7, 2020

Why Choose the Infant Hip Sonography Training Phantom for Your Pediatric Radiology Training Program in India?

The Infant Hip Sonography Training Phantom (Item No. US-13, 3kg product weight, Brand: Kyoto Kagaku, Japan) is described by the manufacturer as the market's only training phantom for hip sonography built on a full-body manikin of a 6-week-old infant — and the best available tool for teaching Graf's method, the internationally recognized classification system for infant hip ultrasound. Production and development were supervised by Univ. Prof. h.c. Reinhard Graf, MD, the method's own originator, alongside the Japanese Society of Orthopedic Ultrasonics.

The full-body infant format facilitates genuine anatomical understanding, and bilateral hips allow for realistic comparative examination — exactly as a sonographer would assess a real infant patient.

Training applications: Setting up and preparing for hip sonography, changing the infant's position during examination, communication and interaction with the infant's guardian, correct transducer use and positioning, recognition of landmarks for hip sonography, visualization of standard, anterior, and posterior planes, and interpretation/morphological classification of the sonogram.

Set includes: Ultrasound infant phantom, talcum powder, storage case, instruction manual.

Specifications: Size W55×D25×H13cm (W21.6×D9.8×H5.1in). Product weight 3kg (6.6lb). Materials: polyurethane elastomer, latex-free.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Kyoto Kagaku products in India with over 30 years of experience supplying medical simulators and pediatric radiology training equipment to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: What is Graf's method? A: The internationally recognized classification system for infant hip ultrasound, used to screen for developmental dysplasia of the hip. This phantom was developed with the method's own originator.

Q: Is this suitable for latex-allergic students or staff? A: Yes — the materials are polyurethane elastomer and explicitly latex-free.

Q: How do I order from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.