



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Cutdown Pad, Light

Item No. ALT00280
Weight
Make: Adam, Rouilly Limited, UK

[Read More](#)

SKU:

Categories:Clinical Skills Trainers



Product Description

A multi-layered pad containing filled veins for practising venous cutdown procedure.

Skills

- Multi-layered pad containing filled veins used for practising the cut-down procedure
- Veins can only be located by palpating surface as they are deep within skin layers
- Tissue incision
- Cannula insertion
- Veins injection both distally and proximally
- Tying off veins. It is recommended veins are tied off each time so that the next student may

repeat the procedure without losing the contents from the veins

Features

- Layered pad simulating the epidermis and subdermal layer
- Three straight 150mm long veins charged with artificial blood
- 10 or 12 cuts can be made on the pad making it economical for training

Important

- Only ALT00223 ARTIFICIAL BLOOD can be used with the cutdown pad

Manufacturer

Limbs & Things