



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com




Model for Gynecological Patient Education - 3B Smart Anatomy

Item No. 1013705 [P53]
Weight 1.091 kg
Dimensions 26 x 19 x 22 cm
Brand 3B Scientific

[Read More](#)

SKU:

Categories: Women's Health Education



Product Description

Trust is one of the core elements of a successful healthcare relationship. It is also one that takes time to build, especially regarding intimate topics. Invest in a long lasting relationship and continued

business and utilize education tools to build a strong rapport with your patients.

This unique gynecological training model is ideal for demonstration purposes and for training realistic insertions of female barrier [contraceptive devices](#), which are placed in the vaginal/cervical area.

These include:

- Femidom
- FemCap (cervical cap)
- Vaginal ring
- Diaphragm

The second important function of this simulator is that it is possible to insert the various types of [pessaries](#) such as ring, bowl or cube which are used in the event of a prolapse or incontinence. It is also possible to demonstrate how to introduce a IUD (coil) and show where it should be located. The cervical cap can be detected with ultrasound and it is also possible to practice intrauterine insemination (IUI).

The transparent vaginal wall allows a view of the inside of the vagina and because the uterus is partially opened up it can be clearly seen how the various devices are inserted. Correct positioning and insertion can be demonstrated and then checked in a teaching scenario, the trainee quickly gains confidence in using the various aids due to a clear representation and simple self-check. This is a very simplified representation of the female pelvis and the female genitalia.