



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



Intubation Head for CPR

Item No. 1019711 [P71/AH]

Weight 1.6kg

[Read More](#)

SKU:

Categories: Airway Management Adult



Product Description

Upgrade CPR*Lilly*PRO™ to include airway management training with the Trucorp® advanced airway head: The head offers the most realistic internal anatomical details for a full range of airway management techniques. It is easy to install on your existing CPR*Lilly*PRO training manikin and can quickly be interchanged when switching between CPR and airway management training scenarios.

Life-like airway for realistic skill training in:

- Double naso-tracheal intubation
- Bag and Mask ventilation techniques
- Full range of supraglottic devices
- Direct laryngoscopy
- Endotracheal tube insertion
- Awake fibre optic examination
- Combi tube insertion

Recommended equipment size for airway management training:

- 7.0-7.5 mm ID for nasal intubation
- 8.0-9.0 mm ID for oral intubation
- Size 3-5 for LMA laryngeal masks
- Similar respective sizes for other supraglottic devices

Manikin can be used with Supraglottic airway devices for airway training. Allergy note: head skin contains latex.