

## **SEM TRAINERS & SYSTEMS**

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile: Email: +91 88495 63724 sem@semtrainers.com Mobile 1: +91 98791 03905 Website: www.semtrainers.com



## Rescue Randy 9000

Item No. 149-9000U

Weight 30.6kg


Make: Nasco Healthcare, USA

**Read More** 

SKU:

Categories: Rescue Training





## **Product Description**

This manikin is designed to be customized by the user for optimal training use. The Rescue Randy 9000 is made from rugged polyethylene parts that allow for the addition of varying amounts of weight to each of the pieces. You may add water, sand, or any other substance as long as you can empty the substance out of the parts. The empty manikin is easily moved and stored when training is completed. This custom-molded construction ensures constant size from one manikin to the next and allows the







utilization of Randy 9000 in every environment imaginable. Each sealed piece is assembled with stainless steel hardware. The range of motion of the joins, including bending at the waist, mimics real life. Randy 9000 is perfect for RIT or Combat Challenge as well as ship-board, heavyweight, horizontal rescues or lightweight packaging exercises like spinal immobilization. Randy 9000 may be used in all weather conditions and the joint design protects rescuer's fingers. Weight with water: 170 lbs Weight with sand: 250 lbs





