



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com

Tracheotomy Care Simulator

Item No. LF01083U
Weight 9.39kg
Dimensions 53 x 28 x 53 cm
Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories:Advanced Trauma Life Support (ATLS)



Product Description

The **Life/form®** Tracheostomy Care Simulator is the most realistic approximation of a patient available today. Careful attention to detail is incorporated into this pseudo-patient. Suctioning techniques, proper cuff inflation, dressing changes, and other techniques can be practiced and critiqued. The simulator replicates an adult male containing oral and nasal passages and all appropriate anatomy - pharynx, epiglottis, trachea, esophagus, stoma, cricoid cartilage, and representative cervical vertebrae. The chest area contains left and right bronchi and the bronchial tree. The esophagus "dead ends" 2" below the trachea opening. The side of the neck and chest area over the bronchial tree allow viewing of the suction and trachea tube. The viewing window at the bottom of the simulator shows the location of a suction catheter in the right bronchus (seen) or left bronchus (unseen). Oral, nasopharyngeal, nasotracheal, and tracheal suctioning can be practiced. Cleansing of the stoma area can be performed just as on an actual patient. Changing dressings and ties is also possible, and the experience is very close to an actual one. Confidence and practice are at your fingertips with this **Life/form®**

Tracheostomy Care Simulator. Comes complete with instruction guide and hard carry case.
Tracheostomy tube is **not** included. Size: 21" x 21" x 11" • Suctioning • Cuff Inflation • Dressing Changes • Visual Landmarks • Visual Verification of Bronchial Location