



**SEM**

## SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

**Mobile :** +91 88495 63724

**Mobile 1:** +91 98791 03905

**Email:** sem@semtrainers.com

**Website:** www.semtrainers.com

## ORTHO

Item No. 1019876  
Weight 0.031 kg  
Dimensions 9.5 x 4.5 x 3.3 cm

[Read More](#)

**SKU:**

**Categories:** 3B Biomechanical Test Blocks



## Product Description

### 3B Biomechanical Test Blocks

- Use to test implants and devices
- Train the evaluation of bone density

- Simulate different sets of characteristics, i. e. bone density and cortical thickness

**Bone density: 10 – 30 PCF (Pounds per Cubic Foot)**

The solid rigid polyurethane foam is available in 4 different densities: 10, 15, 20 and 30 PCF.

**Cortical thickness: 1 – 3 mm**

The cortical shell is simulated by a dense plate on top of each block.

Its density is 40 PCF and three different thicknesses (1, 2 or 3 mm) are available for each foam density option.

**Dimensions:**

CP1: 1 mm cortical thickness: 9.5 x 4.5 x 3.1 cm

CP2: 2 mm cortical thickness: 9.5 x 4.5 x 3.2 cm

CP3: 3 mm cortical thickness: 9.5 x 4.5 x 3.3 cm