



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com

Battery Leg SMART STAT

Item No. 1018576 [W44183]

Weight 4.99kg

Dimensions 78.7 x 20.3 x 20.3 cm

MPN: 101-407NU / PP00407U / SB46528U

[Read More](#)

SKU:

Categories:ALS Adult



Product Description

The SMART STAT Patient Simulator operates by utilizing the battery system in the lower left leg. The battery system is rechargeable and will power SMART STAT for approximately five hours. Charging the battery system take about five hours and having an extra battery leg will allow you to work all day without having to plug into an electrical outlet.