



**SEM**

## SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

**Mobile :** +91 88495 63724

**Mobile 1:** +91 98791 03905

**Email:** sem@semtrainers.com

**Website:** www.semtrainers.com


## Amputated Bleeding Leg

Item No. 1018590 [W44180]  
Weight 2kg  
Dimensions 43.2 x 15.2 x 22.9 cm  
MPN: 101-408U / PP00408U / SB46529U

[Read More](#)

**SKU:**

**Categories:**ALS Adult



## Product Description

Substitute the amputated Bleeding leg for the right leg of the SMART STAT by simply attaching the blood supply line to the knee of a SMART STAT. The instructor controls the bleeding by utilizing the "Venous Blood Flow" on the Pocket PC included with SMART STAT. The SMART STAT holds 500 cc of blood for use with the Amputated Bleeding Leg. A properly applied tourniquet (not included) will arrest the bleeding