



## SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

**Mobile :** +91 88495 63724  
**Mobile 1:** +91 98791 03905

**Email:** sem@semtrainers.com  
**Website:** www.semtrainers.com



## Deluxe Plus CRiSis™ Manikin with CPR Metrix and iPad®

Item No. 1022168  
Weight 28.078kg  
MPN: LF03989[LF03956U]  
Previous Item No. 1005649

[Read More](#)

**SKU:**

**Categories:**ALS Adult

## Product Description

**Life/form® Adult CRiSis™ Manikins** Complete resuscitation system that presents challenges you may face in real world situations, with five major features allowing practice of several simulation scenarios. It is ideal for Advanced Life Support (ALS) courses aimed to teach resuscitation guidelines and skills to healthcare professionals. Modular components allow you to create a manikin to suit your needs. **Airway Management**

- Compatible airway management devices include E.T., E.O.A., E.T.L., L.M.A., E.G.T.A., Combitube®, and King Systems
- Cricoid cartilage allows for practice of Sellick maneuver
- Oral, nasal, and digital intubation
- Palpable carotid pulse (manual)

- Separate left and right lungs for auscultation
- Suctioning capabilities

### **Defibrillation Chest**

- Compatible with all standard brands and types of defibrillators, monitors, and patient simulators
- Internal load box absorbs full strength of every shock — designed to handle up to maximum recommended 360 Joules
- Manual, semiautomatic, and automatic defibrillation
- Monitor manikin like a real patient at 4 ECG sites and 2 defibrillation sites CPR
- Fully articulated head, neck, and jaw
- Jaw thrust
- Palpable and visual landmarks

### **Blood Pressure Arm**

- Auscultatory gap
- 5 Korotkoff phases
- Systolic and diastolic levels, heart rate, and sound volume
- Radial pulse
- Use with an optional external speaker system to broadcast sounds to the entire classroom

### **IV Arm**

- Articulated at the biceps for antecubital and dorsal access
- Bony landmark at shoulder to identify muscle tissue for intramuscular injections
- Realistic flashback from a pressurized system
- Replaceable skin and veins