



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com

Replacement sets for episiotomy trainer (set of 5)



Item No. 1019641 [XP95-002]
Weight 1.84 kg
Dimensions 220 x 130 x 95 cm
Brand 3B Scientific

[Read More](#)

SKU:

Categories:Gynecology



Product Description

Set of 5 silicone inserts for the Episiotomy and Suturing Simulator P95. Each cut can be sutured several times. Muscle and skin layers can be sutured separately. Suturing the following birth injuries

can be practiced:

- Episiotomy, for example:
 - o Median episiotomy
 - o Mediolateral episiotomy
 - o Lateral episiotomy
- Vaginal tears
- Labial tears

Additionally, the following muscle structures can be palpated and sutured:

- M. sphincter ani externus (outer anal sphincter muscle)
- M. transversus perinei superficialis (superficial transverse perineal muscle)
- M. bulbospongiosus (bulbospongiosus muscle)