



## SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

**Mobile :** +91 88495 63724  
**Mobile 1:** +91 98791 03905

**Email:** sem@semtrainers.com  
**Website:** www.semtrainers.com



## Life/form® Child CPR/Airway Management Torso

**Item No.** LF03633U  
**Weight** 3.856kg  
**Make:** Nasco Healthcare, USA

[Read More](#)

**SKU:**

**Categories:**BLS Child

### Product Description

The CPR/Airway Management Torso is ideal for practicing intubation, ventilation, suctioning techniques, and CPR. This torso is designed to be the foundation of the Child CRiSis™ system as arms, legs, and a defibrillation chest skin can be added later to provide IV access, intraosseous infusions, blood pressure monitoring, and ECG recognition and defibrillation training. Includes pump lubricant. Manikin can be used with Supraglottic airway devices for airway training.

Head/torso	Age	Digital intubation	Tongue edema	Laryngospasm	Stomach rise	Suctioning	Manual carotid pulse	Defibrillation	Additional features
------------	-----	--------------------	--------------	--------------	--------------	------------	----------------------	----------------	---------------------

<b>NG Tube Insertion Trainer</b> <i>Item No:</i> <a href="#">1017250</a>	head mounted on a board	1-year old	-	-	-	Yes	-	-	-	trainer for instruction of nasogastric tube insertion and care
<b>Infant Airway Simulator</b> <i>Item No:</i> <a href="#">1018863</a>	full-body	1-year old	-	-	-	-	Yes	-	-	chest cavity containing realistic organs, nasal passage permits placement of NP tube
<b>Infant Airway Trainer</b> <i>Item No:</i> <a href="#">1017954</a>	head mounted on a board	1-year old	-	-	-	Yes	Yes	-	-	accepts 4.0 mm endotracheal tubes
<b>Infant Intubation Head</b> <i>Item No:</i> <a href="#">1017245</a>	head mounted on a board	1-year old	-	-	-	-	Yes	-	-	accepts 2.5 mm endotracheal tubes
<b>Advanced Infant Intubation Head</b> <i>Item No:</i> <a href="#">1017236</a>	head mounted on a board	1-year old	-	Yes	-	Yes	Yes	-	-	Vocal cords are highlighted for easy viewing
<b>Airway Management Trainer</b> <i>Item No:</i> <a href="#">1014201</a>	head mounted on a board	3-year old	-	-	-	-	Yes	-	-	accepts 4.0 mm endotracheal tubes
<b>Advanced Child Airway Management Trainer</b> <i>Item No:</i> <a href="#">1019787</a>	head	5-year old	-	Yes	Yes	-	Yes	Yes	-	can be attached to any of the Life/form® Child Airway Management Trainers or Child CRiSis Manikins

<b>Child Airway Management Trainer Head Item</b> <i>No: <a href="#">1019788</a></i>	head	5-year old	-	-	-	Yes	Yes	Yes	-	can be attached to any of the Life/form® Child Airway Management Trainers or Child CRiSis Manikins
<b>Advanced Child Airway Management Trainer</b> <i>Item No: <a href="#">1005698</a></i>	head mounted on a board	5-year old	Yes	Yes	Yes	Yes	Yes	Yes	-	also allows for E.T., E.O.A., P.T.L., L.M.A., E.G.T.A., Combitube®, and King System insertion
<b>Basic Child CPR/Airway Management Torso</b> <i>Item No: <a href="#">1018865</a></i>	torso	5-year old	-	-	-	-	-	-	-	This torso is the foundation of the Child CRiSis™ system as arms, legs, and a defibrillation chest can be added later.
<b>Advanced Child CPR/Airway Management Torso</b> <i>Item No: <a href="#">1018864</a></i>	torso	5-year old	Yes	Yes	Yes	Yes	Yes	Yes	Yes	also allows for E.T., E.O.A., P.T.L., L.M.A., E.G.T.A., Combitube®, and King System insertion