

## **SEM TRAINERS & SYSTEMS**

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile :	+91 88495 63724	Email:	sem@semtrainers.com
Mobile 1:	+91 98791 03905	Website:	www.semtrainers.com



## Standard Ready-or-Not Tot® - White Male

Item No. Weight Dimensions MPN: Read More 1017931 [W44214] 2.988 kg 56.5 x 29.6 x 16.5 cm WA18306

SKU:

Categories:Neonatal Patient Care



**Product Description** 

The original manikin that simulates a baby's varying needs!

Take the glamour out of teen pregnancy and demonstrate the realities of the real world to your students.

Characteristics of Standard Ready-or-Not Tot®:

• Offers three different tending programs typical of most babies, to simulate normal baby-tending needs along with fussy periods in each program.

• Provides programmed tending events that are unpredictable to the student, but easy for the teacher to monitor, unlike other simulators.

• Cries, coos, and burps in response to appropriate care. Programs are 48 hours in length and run continuously.

- Includes a "demo" option.
- Demands student's attention for periods of 5 minutes to 30 minutes for each tending occurrence. There are 25-27 tending occurrences within each 48-hour time frame (for total of 7 hours of care).
- Provides teacher with complete control of the child care experience.
- Includes a "panic" key for quieting difficulties without discontinuing the entire simulation.
- Alerts teacher to a student's abuse, panic, or tampering with the control box.
- Newborn-size doll represents a baby of a teen mom 5 lbs. 8 oz., 20" long.
- Anatomically correct.

You never have to worry if students neglected the manikin. Ready-or-Not Tot® uses set programs so the teacher always knows when the manikin required care. Student response to the manikin's needs are easily checked using the template provided.

Comes complete with one set of student keys, one set of teacher keys, reproducible parental consent form/permission slip, reproducible student response sheet, teacher correction template, diaper, 9V battery, and instructions.