

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile :	+91 88495 63724	Email:	sem@semtrainers.com
Mobile 1:	+91 98791 03905	Website:	www.semtrainers.com



IV Arm- Child- 5 years old

Item No. Weight Make: Read More LF03612U 3.905kg Nasco Healthcare, USA

SKU:

Categories:Injections and Punctures

Product Description

An exact reproduction of the arm of an 5-year-old for practice of venipuncture and intramuscular injection techniques and procedures used in young children. Features realistic vinyl skin and latex veins so the skin actually rolls as you palpate the vein, and the characteristic "pop" can be felt as the needle penetrates the vein. Realistic flashback occurs with proper insertion. The arm pivots at the upper arm for easy accessibility to antecubital fossa, along the forearm, and at the back of the hand. A bony landmark at the shoulder helps identify muscle tissue for intramuscular injections. For attachment to the left shoulder of the Child CPR/Airway Management Torso (1018865). The arm comes with light skin, two fluid supply bags, 3 cc syringe, 12 cc syringe, needle, one pint Life/form® blood, teaching guide, and carry case.