

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Email: Mobile: +91 88495 63724 sem@semtrainers.com Mobile 1: +91 98791 03905 Website: www.semtrainers.com



IV Arm & Hand Trainer | Injections and Punctures | LF00702U

Item No. LF00702U

Weight 0.8kg

Make: Nasco Healthcare, USA

Read More

SKU:

Categories:Injections and Punctures

Product Description

The economical IV Training Arm and Hand combine realism, fine detail, and lightweight convenience in a single trainer. These self-contained trainers are packaged in a plastic case that can be converted into a work station. Trainers include everything you will need to begin training and practicing IV skills. The Life/form® Portable IV Training Arm and Hand Trainers are constructed of soft material with life-like veins in the skin surface that are visible and palpable. Veins are accessible at the antecubital fossa, along the forearm, the back of the hand, and thumb veins, making it possible to practice venipuncture at any of the common sites. When puncturing through the vinyl skin and veins, the skin will actually roll as you palpate the veins and the characteristic "pop" can be felt as the needle penetrates the vein. These trainers are the ideal way to teach venipuncture techniques including starting IVs and introducing Over the Needle IV catheters, and a realistic flashback will be observed when the needle has been accurately inserted into the vein. These trainers are inexpensive enough so that each student can have their own arm and hand to practice technique and skills required in their courses. Under







normal use, hundreds of injections may be performed. The skins and veins cannot be changed and trainers will eventually need to be replaced. Included with each of the trainers are two IV bags, packet of blood powder and spoon, plastic fluid supply stand, two pinch clamps, 3 cc syringe, 12 cc syringe, 12

gauge needle, a winged infusion set, and a plastic clamshell storage case.

	Phlebotomy	IM site ^S	Subcutaneous injection	Intravenous bolus injection	Sterile technique	Infusion technique ho	Suture, incision and emiodialysis	Pressure squeeze bulb and clenched fist simulation	collapsed veins	Replacement skin and veins
Venatech IV Trainer Item No: 1017966	-	-	-	-	-	-	-	-	-	Skin 1021707
Injection Training Hand Item No: 1017962	-	-	-	-	-	-	-	-	-	on request
Portable I.V. Hand Trainer Item No: 1017958	-	-	-	-	-	-	-	-	-	on request
Portable I.V.Arm Trainer Item No: 1017960	-	-	-	-	-	-	-	-	-	on request
I.V. Training Hand Item No:1005754	-	-	-	-	-	-	-	-	-	on request
I.V. Injection Hand Item No: 1017963	-	-	-	-	-	-	-	-	-	Skin 1005667 Veins 1005665
I.V. Arm & Hand Trainer Item No: 1013746	-	-	-	-	-	-	-	-	-	on request
3B Scientific® I.V. Injection Arm Item No:1021418	-	-	-	,	•	•	-	-	-	Skin 1021426 Veins 1021427







Intravenous Training Arm Item No: 1018755	-	-	-	•	-	•	-	,	•	on request
Injection Arm Item No:1005809	-	in deltoid area	•	•	-	•	-	•	-	Skin 1013384 Veins 1013385
Geriatric IV Arm Item No: 1005777	•	-	-	-	•	•	-	-	-	Skin 1020949 Veins 1021257
SMASH Advance Patient Training Item No: 1019645	•	•	•	•	-	•	•	-	•	Skin 1012333 Veins 1020609





