



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com

Deluxe Child CRiSis™ with Advanced Airway Management

Item No. LF03981U
Weight 11.919kg
Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories:ALS Child



Product Description

Life/form® Child CRiSis™ Manikins This complete resuscitation system is designed for teaching individuals or groups life-saving techniques for children. It provides complete Pediatric Advanced Life Support (PALS) training capabilities and combines all of the skill stations on one manikin, making training more realistic and giving students an appreciation for the smaller scale of the patient. Simulates a 5-year old child, 118 cm tall and 9 kg in weight. **Airway Management**

- Compatible airway management devices include E.T., E.O.A., E.T.L., L.M.A., E.G.T.A., Combitube®, and King Systems
- Inflatable stomach bladder to indicate esophageal intubation
- Oral, nasal, and digital intubation capabilities
- Realistic anatomical features
- Separate left and right lungs for auscultation

- Suctioning capabilities

Intraosseous Infusion

- Palpable landmarks include the patella, tibia, and tibial tuberosity
- Pressurized system allows aspiration of fluid
- Replaceable bones and skin

Blood Pressure Arm

- Auscultatory gap
- Systolic and diastolic levels, heart rate, and sound volume
- Palpable radial pulse
- Speaker in arm reproduces real blood pressure sounds
- 5 Korotkoff sounds

CPR

- Fully articulated head, neck, and jaw thrust
- Visual landmarks and palpable carotid pulse (manual)

IV Arm

- Articulated at the biceps for antecubital and dorsal access
- Realistic flashback from a pressurized system

Defibrillation Chest

- Compatible with all standard brands and types of defibrillators, monitors, and patient simulators
- Internal load box absorbs full strength shock
- Manual, automatic, and semi-automatic defibrillation
- Rhythms are detected at 4 ECG sites and 2 defibrillation sites

Femoral Access

- Palpable arterial pulse (manual)
- Palpable symphysis pubis and anterior superior iliac spine
- Realistic flashback from pressurized venous system

Manikin can be used with Supraglottic airway devices for airway training.

	Basic Child CRiSis™	Basic Child CRiSis™ with Advanced Airway Management	Complete Child CRiSis™	Complete Child CRiSis™ with Advanced Airway Management	Deluxe Child CRiSis™ with ECG	Deluxe Child CRiSis™ with Advanced Airway Management
Item No.	1021992	1021993	1021994	1021995	1005624	1021996
Airway Management	✓	✓	✓	✓	✓	✓
Laryngospasm and tongue edema	-	✓	-	✓	-	✓
Blood pressure arm	-	-	✓	✓	✓	✓
CPR	✓	✓	✓	✓	✓	✓
Defibrillation chest	-	-	✓	✓	✓	✓
Femoral access leg	-	-	✓	✓	✓	✓
Intraosseous infusion	-	-	✓	✓	✓	✓
IV arm	-	-	✓	✓	✓	✓
Interactive ECG simulator	optional 1005666	optional 1005666	optional 1005666	optional 1005666	included	included
Carry case	optional 1005651	optional 1005651	optional 1005651	optional 1005651	optional 1005651	optional 1005651