



## SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

**Mobile :** +91 88495 63724  
**Mobile 1:** +91 98791 03905

**Email:** sem@semtrainers.com  
**Website:** www.semtrainers.com

## Complete Child CRiSis™

**Item No.** LF03616U  
**Weight** 11.919kg  
**Make:** Nasco Healthcare, USA

[Read More](#)

**SKU:**

**Categories:**ALS Child



## Product Description

**Life/form® Child CRiSis™ Manikins** This complete resuscitation system is designed for teaching individuals or groups life-saving techniques for children. It provides complete Pediatric Advanced Life Support (PALS) training capabilities and combines all of the skill stations on one manikin, making training more realistic and giving students an appreciation for the smaller scale of the patient. Simulates a 5-year old child, 118 cm tall and 9 kg in weight. **Airway Management**

- Compatible airway management devices include E.T., E.O.A., E.T.L., L.M.A., E.G.T.A., Combitube®, and King Systems
- Inflatable stomach bladder to indicate esophageal intubation
- Oral, nasal, and digital intubation capabilities
- Realistic anatomical features
- Separate left and right lungs for auscultation

- Suctioning capabilities

### **Intraosseous Infusion**

- Palpable landmarks include the patella, tibia, and tibial tuberosity
- Pressurized system allows aspiration of fluid
- Replaceable bones and skin

### **Blood Pressure Arm**

- Auscultatory gap
- Systolic and diastolic levels, heart rate, and sound volume
- Palpable radial pulse
- Speaker in arm reproduces real blood pressure sounds
- 5 Korotkoff sounds

### **CPR**

- Fully articulated head, neck, and jaw thrust
- Visual landmarks and palpable carotid pulse (manual)

### **IV Arm**

- Articulated at the biceps for antecubital and dorsal access
- Realistic flashback from a pressurized system

### **Defibrillation Chest**

- Compatible with all standard brands and types of defibrillators, monitors, and patient simulators
- Internal load box absorbs full strength shock
- Manual, automatic, and semi-automatic defibrillation
- Rhythms are detected at 4 ECG sites and 2 defibrillation sites

### **Femoral Access**

- Palpable arterial pulse (manual)
- Palpable symphysis pubis and anterior superior iliac spine
- Realistic flashback from pressurized venous system

Manikin can be used with Supraglottic airway devices for airway training.

|                                  | Basic Child<br>CRiSis™              | Basic Child<br>CRiSis™ with<br>Advanced Airway<br>Management | Complete Child<br>CRiSis™           | Complete Child<br>CRiSis™ with<br>Advanced Airway<br>Management | Deluxe Child<br>CRiSis™ with ECG    | Deluxe Child<br>CRiSis™ with<br>Advanced Airway<br>Management |
|----------------------------------|-------------------------------------|--------------------------------------------------------------|-------------------------------------|-----------------------------------------------------------------|-------------------------------------|---------------------------------------------------------------|
| Item No.                         | 1021992                             | 1021993                                                      | 1021994                             | 1021995                                                         | <a href="#">1005624</a>             | 1021996                                                       |
| Airway<br>Management             | ✓                                   | ✓                                                            | ✓                                   | ✓                                                               | ✓                                   | ✓                                                             |
| Laryngospasm<br>and tongue edema | -                                   | ✓                                                            | -                                   | ✓                                                               | -                                   | ✓                                                             |
| Blood pressure<br>arm            | -                                   | -                                                            | ✓                                   | ✓                                                               | ✓                                   | ✓                                                             |
| CPR                              | ✓                                   | ✓                                                            | ✓                                   | ✓                                                               | ✓                                   | ✓                                                             |
| Defibrillation<br>chest          | -                                   | -                                                            | ✓                                   | ✓                                                               | ✓                                   | ✓                                                             |
| Femoral access<br>leg            | -                                   | -                                                            | ✓                                   | ✓                                                               | ✓                                   | ✓                                                             |
| Intraosseous<br>infusion         | -                                   | -                                                            | ✓                                   | ✓                                                               | ✓                                   | ✓                                                             |
| IV arm                           | -                                   | -                                                            | ✓                                   | ✓                                                               | ✓                                   | ✓                                                             |
| Interactive ECG<br>simulator     | optional<br><a href="#">1005666</a> | optional<br><a href="#">1005666</a>                          | optional<br><a href="#">1005666</a> | optional<br><a href="#">1005666</a>                             | included                            | included                                                      |
| Carry case                       | optional<br><a href="#">1005651</a> | optional<br><a href="#">1005651</a>                          | optional<br><a href="#">1005651</a> | optional<br><a href="#">1005651</a>                             | optional<br><a href="#">1005651</a> | optional<br><a href="#">1005651</a>                           |