



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE


SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



Chester Chest™ Standard Arm - Light

Item No. SB52877

Weight

Make: Nasco Healthcare, USA

[Read More](#)

SKU:SB52877

Categories:Intravenous (I.v.) and Arterial



Product Description

The Standard Arm has a dual lumen 5Fr PICC (Peripherally Inserted Central Catheter) exiting the basilic vein from the inner bicep area, which is currently the preferred placement. The basilic vein is slightly raised for easy identification. Also featured is a pre-positioned 20g IV Catheter in the forearm. The dual PICC and 20g IV Catheter are attached to a simulated blood reservoir bag in the arm to permit the practice of "blood" withdrawal and fluid infusion. Cleansing, application of dressings and securement devices can also be demonstrated at all sites. The Standard Arm fits all Chester Chest

models. **Features:**

- Dual lumen 5 FR PICC exiting the basilic vein from the inner biceps area, which is the current preferred placement
- Basilic vein is slightly raised for easy identification
- Pre-positioned IV catheter in the forearm
- Also includes blood bag reservoirs
- 18 in. L