



## SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

**Mobile :** +91 88495 63724  
**Mobile 1:** +91 98791 03905

**Email:** sem@semtrainers.com  
**Website:** www.semtrainers.com

## Life/form® Newborn Nursing Skills and ALS Simulator

Item No. LF01400  
Weight 7-1/2-lb.  
Make: Nasco Healthcare, USA

[Read More](#)

**SKU:**LF01400

**Categories:**Neonatal Patient Care



### Product Description

The Life/form® Newborn Nursing Skills and ALS Simulator is used for meeting objectives of the Neonatal Resuscitation Program's course curriculum. Newborn size: 19 in. L (crown to heel), 7-1/2-lb. weight, 14-1/2 in. head circumference. Includes simulator, replaceable airway (2), bilateral lung, collapsed lung, functional (patent) umbilicus, omphalocele, myelomeningocele, hand skins (2), foot skins (2), vein system for hand, vein system for foot, 500 ml fluid supply bag (2), 3 cc syringe with

needle, 25 gauge infusion butterfly needle, 25 ml syringe, artificial blood powder, lubricant, baby powder, and hard carry case. Airway, Breathing, and Ventilation Anatomical landmarks Auscultation of lung sounds during ventilation Bilateral lung expansion with visually dramatic chest rise ET tube insertion Manually control rate and depth of respiration Positive pressure ventilation Realistic chest compressions Sellick maneuver Unilateral chest movement option (collapsed lung and mainstem intubation) Birth Anomalies Neural tube defect (Myelomeningocele) easily attaches to the back Omphalocele Chest Tube Practice chest tube insertion and care of the chest tube area Cleaning and Diaper Changing Practice cleaning, bathing, and diaper changing Waterproof and completely immersible Gastrointestinal (GI) Accepts a OG/NG tube for practicing tube feeding and suctioning Monitoring Various monitors, sensors, and electrodes may be attached and connected Observation and Measurement Observe and measure each part of the body (head, umbilicus, skin, etc.) Scenarios Build your own high-risk delivery, emergency, transport, and NICU scenarios Skin Delicate appearance and texture of real skin Soft, durable lifelike skin Suction Practice suctioning nose, oral cavity, and stomach Umbilicus Can be attached and detached easily Made of soft, realistic material Practice blood withdrawal or infuse fluids Patent umbilicus with venous and arterial access Venous Access Accepts IV catheter or butterfly infusion needle IV hand and foot with realistic flashback from pressurized system Soft, lifelike skins with small lumen tubing for realistic practice 5-year warranty.