

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Email: Mobile: +91 88495 63724 sem@semtrainers.com Mobile 1: +91 98791 03905 Website: www.semtrainers.com



Chester Chest™ with Peripheral Port Access Arm

Item No. SB18636

Weight

Brand Nasco Healthcare

Read More

SKU:

Categories:Intravenous (I.v.) and Arterial

Product Description

Lifelike model of a human torso Can be used upright or in a supine position Also has recessed area posterior to PICC insertion site for optional peripheral port (peripheral port not included) Base of recessed area is made of soft material that permits port to "float" when accessed When placed over port, included tissue flap allows for a realistic feel when palpating and accessing External Central Catheter Pre-positioned, surgically placed in the right chest Subcutaneous tunnel is visible up to and just over the clavicle Distinguishable Dacron® cuff Distal catheter end is pre-attached to the blood reservoir bag Implanted Vascular Access Device Located in the left chest under a specially formulated flap designed to simulate real human tissue Situated over a rigid underlying surface with molded ribs and a recessed area for the interchangeable inserts 3 difficult accessing inserts placed either under or over the port to simulate palpating and accessing a port with one of the following types of placements: normal, "tipping," "wandering," or "deeply placed" Peripheral Placed Central Catheter Located in the rotated, extended, detachable right arm Dual lumen 5FR PICC exits the slightly raised, basilic vein on the right inner biceps Pre-positioned 20G IV catheter in the right forearm Catheters are pre-attached to







a blood reservoir bag Complete blood draws, heparinization, and fluid infusion while attached or unattached to the torso Additional Capabilities: Subclavian catheter Triple lumen catheter Jugular catheter Includes: Life-size adult torso with a base Tunneled central catheter 9.6 FR Outer tissue pad for left chest Real port (IVAD) 3 difficult accessing inserts Arm with greater degree of rotation and extension than previous arms Dual lumen PICC 5FR Includes user manual 20 G IV catheter Blood reservoir bags Talc (cornstarch) 20-1/2 in. x 15-1/2 in. x 5-1/4 in. 1-year limited warranty