



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com




Bone Fracture Pediatric Phantom "PBU-70B"

Item No. PH-2D
Weight 20 kg
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories: Radiology Phantom



Product Description

Improve skills in detecting bone fractures in children and cultivate awareness of child maltreatment.

Features

Training skills / Applications

1. Training in pediatric radiography can be enriched with clear and subtle bone fractures 2. Typical fractures resulting from child abuse are also included

Plain X-ray CT Basic patient positioning

Case / Pathology

Set includes

Height (approx.)

Packing size (approx.)

Weight (approx.)

Packing weight (approx.)

Materials

Replacement parts

Optional parts

Related products

Update

Signs of callus / Supracondylar humerus fracture / Back, scapula and rib fractures/ Skull fractures
1 pediatric whole body phantom / 1 head supporter / 1 hand fixture belt / 1 screw driver / 1 set of sample X-ray data(DVD) / manual
110 cm / 43.3 in
W86 x D60 x H32cm / W33.9 x D23.6 x H12.6in
20kg / 44lb
25kg/ 55lb (41350-500), 44kg /97lb (41350-600)
Soft tissue: urethane based resin (density: 1.06) / Synthetic bone: epoxy resin (density: 1.31) / Skull: epoxy resin (density:1.11) *Phantom has no metal parts or liquid structure
41350-500-01 head for PH-2D / 41350-500-02 left hand for PH-2D / 41350-500-03 left upper arm for PH-2D / 41350-500-04 torso for PH-2D / 41350-500-05 left thigh for PH-2D / 41350-500-06 left lower thigh-foot for PH-2D
41363-080 storage case
41350-300 PH-2C Pediatric Whole Body Phantom "PBU-70" 41912-100 PH-50 Newborn Whole Body Phantom "PBU-80" 41337-300 Pediatric Chest Phantom 41350-000-11 Fractured Hand / Forarm Phantom
July 28, 2020