



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com




Item No. PH-61
Weight
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories: Radiology Phantom



Product Description

Sectional phantoms allow for imaging of individual anatomy as needed

Features

1. Opaque and transparent types for a diverse training possibility
Opaque: advanced version with close-to-reality challenges in imaging
Transparent: visible bones facilitate understanding in keys for positioning
2. Movable joints of the knee and the elbow for realistic positioning

Set includes

Materials

List of PH-61 series

Update

(For each) 1 phantom / 1 set of sample X-ray data (DVD) / manual
Soft tissue: urethane based resin / Synthetic bone: epoxy resin (density: 1.31) / Skull: epoxy resin (density: 1.2) / *Phantom has no metal parts or liquid structure
41926-010 Head(Transparent) 41926-060 Thorax(Opaque) 41926-070 Thorax(Transparent) 41926-080 Pelvic(Opaque) 41926-140 Right Elbow(Opaque) 41926-150 Right Elbow(Transparent) 41926-020 Right Hand(Opaque) 41926-030 Right Hand(Transparent) 41926-040 Left Hand(Opaque) 41926-050 Left Hand(Transparent) 41926-180 Right Knee(Opaque) 41926-190 Right Knee(Transparent) 41926-100 Right Foot(Opaque) 41926-110 Right Foot(Transparent) 41926-120 Left Foot(Opaque) 41926-130 Left Foot(Transparent)
July 27, 2020