



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com


Subsolid Nodules Phantom

Item No. PH-58
Weight
Brand KYOTO KAGAKU, JAPAN
[Read More](#)

SKU:

Categories: Radiology Phantom

Mixed GGO with single concentric solid field (No.1-7)




Item No.	GGO field		Solid field		Type
	Diameter	HU	Diameter	HU	
1	1.5 cm/0.59 in	-650	0.5 cm/0.20 in	-50	Concentric
2				0	
3				50	
4	2.0 cm/0.79 in		0.3 cm/0.12 in	0	
5			0.5 cm/0.20 in		
6			0.7 cm/0.28 in		
7			0.9 cm/0.35 in		

Mixed GGO with single eccentric solid field (No.8-10)



Item No.	GGO field		Solid field		Type
	Diameter	HU	Diameter	HU	
8	1.5 cm/0.59 in	-650	0.5 cm/0.20 in	-50	Eccentric
9			0	0	
10			50	50	

Mixed GGO with double eccentric solid fields (No.11, 12)



Item No.	GGO field		Solid field		Type
	Diameter	HU	Diameter	HU	
11	2.0 cm/0.79 in	-650	0.3 cm/0.12 in	0	Eccentric
			0.5 cm/0.20 in		
			0.7 cm/0.28 in		
12	2.0 cm/0.79 in	-650	0.5 cm/0.20 in	0	Eccentric
			0.7 cm/0.28 in		
			0.9 cm/0.35 in		

Pure GGO (No. a-h)




Item No.	GGO field	
	Diameter	HU
a	1.5 cm/0.59 in	-750
b		-650
c		-550
d		-450
e		-350
f		-250
g		-150
h		-50

3D GGO



Item No.	GGO field	
	Diameter	HU
3D-GGO	1.5 x 1.5 cm (0.59 x 0.59 in)	-590



Product Description

Both mixed and pure GGO are provided in a variety of sizes and HU numbers.

Features

Case / Pathology

Related products

Update

1. Subsolid Nodules Phantom is a set of simulated lesions designed for study and training in Grand-Glass Opacity (GG0) detection and interpretation. 2. The set also includes 3-D GG0 modeled on clinical CT data. The simulated lesions can be attached to the pulmonary vessels of the Chest Phantom N1 "LUNGMAN" or in the CT Lung Phantom.
41923-000: 7 x Mixed concentric GG0 41923-100: 8 x Mixed eccentric GG0 41923-200: 2 x Mixed eccentric GG0 (double solid fields) 41923-300: 8 x pure GG0 41923-400: 1 x 3D GG0
PH-1 Multipurpose Chest phantom N1 "LUNGMAN"
July 25, 2020