



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Multi Energy CT Quality Assurance Phantom (TR-J)

Item No. PH-75A
Weight 15.5 kg
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories: Radiology Phantom



Product Description

Water Equivalent Material, various inserts and empty bottles enable to verify the appropriate Multi-Energy CT settings

Features

Training skills / Applications

Set includes

Size (approx.)

Packing size (approx.)

Weight (approx.)

Packing weight (approx.)

Materials

Optional parts

Related products

Production ? Development Supervision

Update

1.Phantom using innovative Water Equivalent Material / 2.Empty bottles enable to put various items and check how they / 3.react to Multi Energy/ 4.Save time and efforts to produce custom-made water phantoms / 5.Two different sizes of body. (TR-I, TR-J)
ME-CT image analysis protocol / Metal artifact reduction / Reduction of contrast media
1 truck phantom / 1 internal cylindrical phantom / 9 filling rods for dosimeter holes / 8 WEM inserts / 3 iodine inserts (4, 8, 12mgI/mL) / 1 titanium insert / 1 soft tissue insert / 1 water container insert/ 20 empty sample bottle / 8 spacers for sample bottle / manual
[TR-J] 30×D18×H20cm
W44 x D36 x H42cm / W17.3 x D14.2 x H16.5in
TR-J : 15.5kg
W30.0 x D20.0 x H18.0cm/W11.8 x D7.9 x H7.1 in
Polyurethane
Angle adjustment holder
Multi Energy CT Quality Assurance Phantom (TR-I))
Katsuhiro Ichikawa, Ph.D., / Institute of Medical, / Pharmaceutical and Health Sciences, / Kanazawa University
October 25, 2024