



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905




Email: sem@semtrainers.com
Website: www.semtrainers.com

Multi Slice CT Phantom MHT

Item No. PH-9
Weight
Brand KYOTO KAGAKU, JAPAN
[Read More](#)

SKU:

Categories: Radiology Phantom



Product Description

The phantom can be used for features of CT evaluation such as high and low contrast resolutions, feed direction and CTDI

Features

Training skills / Applications

Set includes

Materials

Optional parts

Update

1. Non-aqueous/Easy Set-up enables liquid-free evaluation session 2. The phantom is designed to allow evaluation in volume scanning
CT
1 low contrast phantom / 1 high contrast phantom / 1 elliptical absorber / 1 low contrast phantom with CTDI / 1 micro disc phantom / 1 Angle adjustment holder (table-top type) / storage case / manual
Acrylic resin, polyurethane
41334-110 sliding phantom holder
August 10, 2020