



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com


Ladder Phantom

Item No. PH-9-2
Weight
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories: Radiology Phantom



Product Description

Evaluation of spatial resolution if simulated contrast enhanced vessels in CT

Features

Training skills / Applications

1. On each plate phantom of 5 mm thickness, five slits of 5 mm length are made to represent vessels 2. Nine variations of vessel width are prepared

CT

Set includes

Materials

Optional parts

Update

1 Outer phantom / 9 ladder phantoms / 1 storage case
Measurement region: epoxy resin, hydroxyapatite / Base: acrylic resin
41919-010 Angle adjustment holder (table-top type)
August 10, 2020