



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com


CT-DI Phantom (Head and Body Phantom)

Item No. PH-59
Weight
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories: Radiology Phantom



Product Description

A set of phantoms for CTDI-100, conforming to requirements described in 21 CFR 1020.33, IEC 61223-3-5: 2004, and IEC 61223-2-6: 2006 as consistency test

Features

1. Represent adult head and body as well as pediatric body
2. This phantom can be used for initial and follow-up QA tests described in JIS Z 4752-3-5:2008 and JIS Z 4752-2-6:2012

Training skills / Applications
Set includes
Size (approx.)
Weight (approx.)
Materials
Related products
Update

Computed Tomography Dose Index (CTDI) Dose profile
1 head phantom / 1 body phantom / 10 filling rods / manual
Body phantom: 32 dia. x 15 cm/12.6 dia. x 5.9 inch Head phantom: 16 dia. x 15 cm/6.3 dia. x 5.8 inch
Body phantom:15kg/ 33lb Head phantom:4kg/ 8.7lb
Acrylic resin
PH-54 CT QA Phantom JCT II
August 10, 2020